

FOLDING STAIR AND EVACUATION CHAIRS

- **34068 STAIR STRETCHER - 131x50x155 cm**
 - **34069 ELECTRIC STAIR STRETCHER - 123x50x160 cm**
- Reliable and safe stretchers with ergonomic design, aluminium alloy frame, easy to handle. Used at home for old people, school security, hospitals, hotels etc. for emergency evacuation. Both climb up and down stairs smoothly and steadily by rubber track. They can move on the ground as manual wheelchairs.

Strong grip track prevents from gliding or falling down when it stops on the stairs. Foldable stretchers, easy to carry and store.



Maximum load: 100 kg. Equipped with 4 wear-resistant rubber wheels, foldable handles, telescopic handle, detachable pillow, PU sponge mattress and safety belts. Height adjustable backrest (34068 only). The electric model (34069 only) has an electric control system that provides:

- circuit inside the pipe, waterproof and noiseless
- LCD display with LED night light
- compliance with EU safety standard



AUTOCLAV. RESUSCITATORS



CE 0197

SILICONE AUTOCLAVABLE RESUSCITATORS - WITH HANDLE
Transparent resuscitators made of medical grade silicone and fully autoclavable at 121°C. Integrated handstrap for user comfort and complete control. Made with no natural rubber or latex. Multilingual manual: GB, FR, IT, ES, PT, SE, RO, GR. See accessories at page 370.

GIMA code	RESUSCITATORS			RESUSCITATORS KIT	
	34260	34261	34262	34263	34264
Type	Adult	Child	Infant	Adult	Child
Mask	n° 5	n° 3	n° 1	n° 5	n° 3
Pressure relief valve H ₂ O	60 cm	40 cm	40 cm	60 cm	40 cm
Volume ml	1,600	550	300	1,600	550
Reservoir bag ml	-	-	-	2,500	2,500
2 m oxygen tube	-	-	-	yes	yes
Reservoir valve	-	-	-	yes	yes

VALUE EVACUATION CHAIR

• **34073 VALUE EVACUATION CHAIR - 4 wheels**
Stair stretcher mainly used to transfer patients up and down stairs or other narrow areas. Structure made of high-strength aluminium while surface is made of leather fabric which can be removed by zipper on the back. Equipped with 4 wheels, four foldable handles on the back and belts for patients safety during transport.

TECHNICAL SPECIFICATIONS

Max. load: 159 kg
Weight: 8 kg
Size: backrest 91x50 cm
seat 49x50 cm
Folded size: 92 x 52 x 16 cm



LED DERMATOSCOPE 10X

Supplied in a rigid case with contact plate and multilingual manual: GB, FR, IT, ES, PT.



- **31186 GIMA POCKET LED DERMATOSCOPE 10X**
Modern, compact and lightweight LED pocket dermatoscope with 10X magnification for evaluation of cutaneous lesions. Precision 10X optical system with achromatic lenses provides distortion free sharp image. Graduated skin contact plate to measure pigmented skin lesions offers accurate diagnoses. Focal range: ± 5D. Focal length: 25 mm. It can be powered by two 1.5 AA alkaline batteries or rechargeable battery (both not included)
- **31586 CHARGING STATION**
- **31587 RECHARGEABLE BATTERY Li-ion 2.5 V**