

GERATHERM® CLASSIC AND CLASSIC+ THERMOMETERS

Mercury-free glass thermometers developed by Geratherm® Medical AG, containing a patented liquid alloy consisting of gallium, indium and tin-called Galinstan. Galinstan is completely non-toxic and environment-friendly. They are ideal for determining body temperature of people of all ages, suitable for rectal, oral or axillary measurement. Each thermometer is calibrated to guarantee measurement accuracy.

- measuring range: 35.5°C to 42°C
- sustainable compared to digital thermometers
- easy to disinfect thanks to the closed glass body, ideal for allergy sufferers
- reset tool included

"Easy flip" shaking aid

Thermometers include a reset tool "easy flip" which makes it easier to lead back the measuring liquid of the thermometers by considerably less effort when shaking down.

Storage case with enlargement function (25532-33 only)

Special storage case with magnifying glass to help people read results more clearly.



25532

Supplied with individual box protective case, reset tool, storage case with magnifying lens (25532-33 only).and multilingual manual. Thermometers are sold in multilingual display of 24 pieces.



25529, 25530



25532, 25533

MADE IN GERMANY



GIMA code	CLASSIC THERMOMETERS display of 24 boxes
25529	Classic therm - GB, FR, IT, ES, PT, DE, NL, NO, DK, FI, SE
25530	Classic therm - GB, PL, LT, LV, EE, CZ, SK, SI, HU, HR, RS, RO
GIMA code	CLASSIC+ THERMOMETERS display of 24 boxes
25532	Classic+ therm - GB, FR, IT, ES, PT, DE, NL, NO, DK, FI, SE
25533	Classic+ therm - GB, PL, LT, LV, EE, CZ, SK, SI, HU, HR, RS, RO

INSECT BITE HELPER

• 24394 INSECT BITE HELPER

Compact and portable device suitable for quick relieving itching, swelling or inflammation caused by insect bites or stings. During application, heat is applied to the treatment area; this can have a positive effect on the healing process.



- dual-switch over-temperature protection design, prevents skin burns due to high temperature
 - two modes with three temperature grades each. The set grade can be automatically saved
 - indicator lights for temperature level
 - low-battery indication
 - skin friendly and chemical free
- Supplied with USB cable and user manual in GB, FR, IT, ES, DE, PL.



RIESTER SUPERVU AND SUPERVU HiRES GALILEAN BINOCULAR LOUPES

SUPERVU GALILEAN BINOCULAR LOUPES

They offer 2.0, 2.5 or 3.0 X magnification which allows greater distance between the user and the patient. They provide a carefully calculated blend of comfort, clarity and field of view.

- excellent optical features with respect to the field view ($\pm 10\%$)
- movement whilst maintaining focal clarity is possible due to superior depth of field.

The lightweight aluminium bar with individual pupillary adjustment provides the mounted loupes with stability and optimal focal alignment, eliminating headache and eyestrain

- ideal for ocular and dental examinations as well as plastic surgery and veterinary professionals.



30965+
30950-57

99 g (including frame).

SUPERVU HiRES GALILEAN BINOCULAR LOUPES

They provide a high-resolution loupe, which far surpasses any previous standards. This 5-element loupe offers triple magnification with a 15% increased field of view, a bright image and a visually flat field from edge to edge.

- excellent optical features with respect to the field view ($\pm 10\%$)
- 3.0 X magnification
- excellent distinction, stunning performance as well as great looks
- ideal for most fields of surgery, surgical assistants as well as ophthalmologists



30966+
30960-62



GIMA code	Riester code	SUPERVU GALILEAN BINOCULAR LOUPES	Magnif.	Working distance	Field of view Ø
30950	3960-620	Binocular loupe	2.0 X	340 mm	114 mm
30952	3960-625	Binocular loupe	2.5 X	340 mm	81 mm
30953	3961-625	Binocular loupe	2.5 X	420 mm	100 mm
30956	3960-630	Binocular loupe	3.0 X	340 mm	59 mm
30957	3961-630	Binocular loupe	3.0 X	420 mm	71 mm
FRAME AND BAND FOR BINOCULAR LOUPES - optional					
30965	12770	Spectacle frame, black without loupe			
30966	12771	Brow band without loupe			

GIMA code	Riester code	SUPERVU HiRES GALILEAN BINOCULAR LOUPES	Magnif.	Working distance	Field of view Ø
30960	3970-630	Binocular loupe	3.0 X	340 mm	68 mm
30961	3971-630	Binocular loupe	3.0 X	420 mm	82 mm
30962	3972-630	Binocular loupe	3.0 X	460 mm	88 mm
FRAME AND BAND FOR BINOCULAR LOUPES - optional					
30965	12770	Spectacle frame, black without loupe			
30966	12771	Brow band without loupe			