













HEINE ORIGINAL BULBS

HEINE BULBS	HEINE ref.	GIMA code	Compatible with following HEINE products
 Xenon 2.5 V	035	31769	Transilluminator Finoff alpha [®] F.O. and F.O. SP Laryngoscopes - Sanalon handle Glaucotest
 Xenon 2.5 V	037	31771	K100 and Beta [®] 100, Beta [®] 100 vet otoscopes Laryngeal mirrors F.O. mini 2000 otoscope Alpha [®] and Alpha F.O. otoscopes Mini/mini 3000/Alpha/Beta tongue holders
 Xenon 3.5 V	049	31787	F.O. Alpha [®] , Beta 100, K 100, Beta 100 vet otoscopes; Alpha [®] , Beta [®] tongue holders
 Miroflex ophthalmoscope	038	31773	Miroflex ophthalmoscope HFR2 spot retinoscope Delta [®] 10 dermatoscope
 Fibrilux mini otoscope	041	31774	Fibrilux mini otoscope Mini 2000 combilamp Mini 1000, 2000 cliplamp
 Alpha ophthalmoscope	042	31780	Alpha [®] ophthalmoscope Mini 2000 ophthalmoscope Alpha [®] Focalux, Mini 2000 Focalux Mini Miroflex ophthalmoscope
 Minilux otoscope	056	31770	Minilux otoscope Mini 2000 otoscope
 Beta ophthalmoscope	069	31781	Beta [®] 200, 200 M2, 200S ophthalmoscopes - XHL 2.5 V
 Beta ophthalmoscope	077	31772	Beta [®] 200, 200 VET, 400; K 180 otoscopes Lambda 100 retinometer Delta [®] 10 plus dermatoscope
 Beta ophthalmoscope	070	31783	Beta 200, 200 M2, 200S ophthalmoscopes - XHL 3.5V
 Beta ophthalmoscope	078	31784	Beta 200, 200 Vet, K180 otoscopes lambda 100 retinometer otoscopes - XHL 3.5V
 K180 ophthalmoscope	084	31782	K180 [®] ophthalmoscope - 2.5 V
 Mini 3000 F.O. otoscope	105	31776	Mini 3000 F.O. otoscope
 Mini 3000 ophthalmoscope	106	31777	Mini 3000 ophthalmoscope
 Mini 3000 combilamp	107	31778	Mini 3000 combilamp and cliplamp - XHL 2.5 V
 Mini 3000 dermatoscope	109	31779	Mini 3000 dermatoscope
 Mini 3000 otoscope	110	31775	Mini 3000 otoscope
 Alpha [®] , Beta 200 retinoscopes	089	31786	Alpha [®] , Beta 200 retinoscopes
 LED 2.5 V	130	31788	Beta 200 LED/halogen retinoscopes

WELCH ALLYN BULBS

WELCH ALLYN BULBS	W. A. ref.	GIMA code	Compatible with following WELCH ALLYN products
 Anoscope	00200-U	32120	Anoscopes: 38800, 380, 385, 395 Illuminators: 41000, 26030 Old products: 235, 77900, 26000, 20100, 21600, 42700, 40510
 Ophthalmoscope	03000-U	32123	Ophthalmoscope: 11710, strabismoscope: 12400, episcope: 47300 Old ophthalmoscope: 11600, 11605, 11610 Retinoscope: 18000
 Handle	03100-U	32124	Otoscopes: 25020, 21700, 20200 Handle: 73500, holder: 28100 Illuminators: 27000, 27050, 26530, 41100, 43300 Old otoscopes: 25000, 25200, 20000
 Ophthalmoscope	03300-U	32125	Ophthalmoscope: 11511-11500
 Otoscope	03400-U	32126	Otoscopes: 24011-24020-24000-24031-21110/1 Handle: 73550, Illuminators: 27200-27250-41110
 Pan optic	03800-U	32127	Pan optic: 11810-11820-11800
 Exam light	04100-U	32128	Exam light: 48400-48410
 Exam light	04200-U	32129	Exam light: 48600, 48610-48625-48635, 48700
 Laryngo	04700-U	32131	Laryngo: 680, 690, 692 (size 1-2) Old laryngo: 634, 664, 674, 684, 694 (check size)
 Old laryngo	04800-U	32132	Laryngo: 680 (size 2-4), 690, 692 (size 3-4) Old laryngo: 634, 640, 664, 674, 684, 694 (check size)
 Ophthalmoscope	04900-U	32133	Ophthalmoscope: 11720, 11730, 11735, 11620, 11630
 Laryngo handles	06000-U	32134	Laryngo handles: 60813, 60814, 60713, 60815, 60803, 60804
 Macroview otoscope	06500-U	32135	Macroview otoscope: 23810, 23820
 Handle	07800-U	32136	Handle: 73211, 73210, 78000, 73200

RIESTER ORIGINAL BULBS

RIESTER BULBS	RIESTER ref.	GIMA code	Bulb	Compatible with following RIESTER products
 2.7 V	10421	31850	Vacuum 2.7 V	Otoscope uni [®] /speculight
 2.5 V	10424	31851	XL 2.5 V	May-ophthalmoscope uni [®]
 3.5 V xenon	10487	31825	XL 3.5 V xenon	Otoscope ri-scope [®] L1
 2.7 V	10488	31852	Vacuum 2.7 V	Pen-scope [®] and e-scope [®] otoscope
 2.5 V XL	10489	31560	XL 2.5 V	Otoscope pen-scope [®] /ri-scope [®] L1 and e-scope [®]
 2.5 V XL	10590	31853	XL 2.5 V	Otoscope uni [®]
 2.5 V XL	10600	31854	XL 2.5 V	Otoscope ri-mini [®] , ri-scope [®] L2/L3 and e-scope [®]
 2.5 V HL	10605	31562	XL 2.5 V	Ophthalmoscope ri-mini/ri-scope [®] L1, L2, L3, e-scope [®] , ri-derma [®]
 3.5 V XL	10607	31855	XL 3.5 V	Ri-scope [®] otoscope L2, L3
 3.5 V XL	10608	31826	XL 3.5 V	Ri-scope [®] ophthalmoscope L1, L2, L3
 LED 3.7 V	14041	31561	LED 3.7 V	E-scope F.O otoscope
 LED 3.7 V	14051	31563	LED 3.7 V	E-scope ophthalmoscopes



COMPATIBLE GERMAN BULBS

GIMA BULBS	GIMA code	Compatible with the following products
LARYNGOSCOPE BULBS - MADE IN GERMANY		
 Conventional laryngoscope SMALL bulb 2.7 V	34324	Gima, Riester, And, Kawe, Timesco Truphatek and most of conventional laryngoscopes in the market (generally used for blades N ^o 00, 0, 1)
 Conventional laryngoscope LARGE bulb 2.7 V	34325	Gima, Riester, And, Kawe, Truphatek, Timesco and most of conventional laryngoscopes in the market (generally used for blades N ^o 2, 3, 4)
 Xenon halogen bulb for F.O. Laryngoscope and otoscope	34485	Gima, And, Kawe and most of F.O. laryngoscopes in the market. Compatible for products that use Heine Bulb # 035
OTOSCOPE-OPHTHALMOSCOPE BULBS - MADE IN FRANCE-GERMANY		
 Xenon halogen otoscope bulb	31795	Heine mini 3000 otoscopes
 Xenon halogen F.O. otoscope bulb	31796	Heine mini 3000 F.O. otoscopes
 Otoscope bulb 2.5 V	31449	Parker otoscopes and some otoscopes of Riester, Timesco and Kawe
 Ophthalmoscope bulb 2.5 V	31446	Parker ophthalmoscope Parker diagnostic sets
 Xenon halogen bulb for otoscope 2.5 V	31478	Gimalux otoscopes, Kawe and Riester mini otoscopes Compatible for products that use Heine bulb #037
 Halogen bulb for ophthalmoscope and dermatoscope 2.5 V	31479	Parker halogen ophthalmoscope
 Xenon halogen ophthalmoscope bulbs 3.5 V	31428	Gima Xenon - halogen diagnostic sets - ophthalmoscope
OTO-OPHTHALMO-DERMATOSCOPE BULBS - OTHER COUNTRIES		
 Xenon halogen otoscope bulbs 3.5 V	31427	Gima Xenon - halogen diagnostic sets - otoscope