

GLUCOSE MONITORS, mmol/L or mg/dL UNITS - 24 LANGUAGES

ADVANCED FEATURES FOR DIABETES MANAGEMENT:

- No Coding** no-coding technology
- HCT** 10-70% HCT range allows neonatal and anemic blood samples
- GDH-FAD** enzyme allows capillary, venous, arterial blood samples. No interference from maltose, galactose...
- Large display screen** for easy reading
- Strip ejector** avoid cross-contamination
- 1,000** test memory with date and time
- Auto sample identification** (blood/control solution)
- 7, 14, 30, 60 and 90-day averages calculation**
- Ketone & Hyper/Hypo** Ketone & Hyper/Hypo warnings
- Meal-marker** function to monitor results before and after meals
- Automatic detection** of insufficient sample



Powered by ACON

Available in mmol/L or mg/dL

5 YEARS WARRANTY

UK REP

Only 24103, 24108, 24109, 24110, 24112, 24114, 24118, 24119, 24120, 24121

CH REP

Only 24102, 24119

GIMA GLUCOSE MONITORS

Accurate and reliable glucose monitors for self testing and professional use. Extensive clinical studies demonstrate the accuracy of Gima glucose monitoring system, powered by Acon, with fresh capillary, venous, neonatal blood samples. Comply with ISO15197:2013/EN ISO15197:2015.

Available as meter with lancet device or complete kit.

Kit includes 10 strips, control solution, 10 sterile lancets, lancet device, clear cap (for testing on forearm or palm), carrying case, user manual.

Bluetooth version (24114 only):

Fully integrated Bluetooth allows for data transmission to data management platforms.

Available in 24 languages.

Please refer to the table below to select the following languages: GB, FR, IT, ES, PT, DE, NL, DK, SE, FI, LT, LV, EE, PL, CZ, HU, SI, HR, RO, BG, GR, Arabic.

GIMA code	GIMA GLUCOSE MONITOR AND ACCESSORIES	Box and manual languages
GIMA GLUCOSE MONITOR mmol/L		
24102	Glucose monitor kit	FR, IT, DE, LV, EE, FI
24103	Glucose monitor kit	GB, NL, SI, LT, HU, BG
24104	Glucose monitor kit	RO, DK, SE, CZ, HR, Arabic
GIMA GLUCOSE MONITOR mg/dL		
24108	Glucose monitor with lancet device	GB, IT, SE, FI
24109	Glucose monitor with lancet device	GB, FR, IT, ES, PT, DE
24110	Glucose monitor kit	GB, FR, IT, ES, PT, DE
24111	Glucose monitor kit	IT, DE, GR, Arabic
24112	Glucose monitor kit	GB, IT, PL, HU, RO, GR
24114	Glucose monitor kit with Bluetooth	GB, FR, IT, ES, PT, DE
ACCESSORIES AND CONSUMABLES FOR GIMA GLUCOSE MONITOR		
24118	Glucose Strips - box of 25 strips*	GB, FR, IT, ES, SI, BG
24119	Glucose Strips - box of 50 strips*	GB, FR, IT, ES, DE, LT
24122	Glucose Strips - box of 50 strips*	PT, NL, SE, CZ, RO, GR
24123	Glucose Strips - box of 50 strips*	PL, DK, FI, HU, HR, Arabic
24120	Glucose Strips - box of 100 strips*	GB, FR, IT, ES, EE, LV
24121	Control Solution	GB, FR, IT, ES
23919	Auto lancet device	
23916	Lancets 28G for 23919 - sterile - box of 100 pcs	
23918	Lancets 30G for 23919 - sterile - box of 100 pcs	
24125	USB cable Glucose Monitor - PC	

*Maximum shelf life: closed vial 24 months, open vial 18 months.



24118



24114

TECHNICAL SPECIFICATIONS	
Technology:	Biosensor/Electrochemical, Glucose dehydrogenase (GDH-FAD)
Calibration method:	No coding
Result calibration:	Plasma-equivalent
Test time:	5 seconds
Sample size:	0.6µL
Sample type:	Fresh whole blood (capillary, venous, neonatal)
Test sites:	Fingertip, forearm, palm
Hematocrit range:	10-70%
Memory storage:	1,000 results with date and time
Test averaging:	7, 14, 30, 60, 90-day averages
Data transfer:	USB cable - optional
Control solution:	3 levels
Automatic shutoff:	2 minutes after last action
Battery:	2x CR 2032 3.0V coin cell batteries
Battery life:	3,000 measurements
Other functions:	Automatic detection of insufficient sample 2nd sample application allowed Sample perturbation detection