


FIAB ECG UNIVERSAL ELECTRODES



All with silver chloride sensor and S/S snap

UNIVERSAL DISPOSABLE ECG ELECTRODES - latex free									
GIMA code	FIAB code	Material	Size	Gel	Shape	Connection	Patient	Application	Minimum order
33311	F9046	PE-Foam	28x44 mm	Solid	Rectang.	Snap	Adult, Pediatric	rest, stress, ECG holter	box of 2,000 (20 bags of 100)
33312	F9046 RM	PE-Foam	28x44 mm	Solid	Rectang.	Carbon Snap	Adult, Pediatric	rest, stress, ECG holter, MRI	box of 2,000 (20 bags of 100)
33313	F9089P-100	Foam	Ø 32x36 mm	Solid	Oval	Snap	Pediatric	rest, stress, ECG holter	box of 2,000 (20 bags of 100)
33314	F9089-100	Foam	Ø 36x40 mm	Solid	Oval	Snap	Adult	rest, stress, ECG holter	box of 2,000 (20 bags of 100)
33370	F9060P	PE-Foam	Ø 40 mm	Solid	Round	Snap	Pediatric	rest, long term, stress, ECG holter	bag of 50
33371	F9060	PE-Foam	Ø 48x50 mm	Solid	Oval	Snap	Adult	rest, long term, stress, ECG holter	bag of 50
33344	F9067	PE-Foam	Ø 48x50 mm	Liquid	Oval	Snap	Adult	rest, stress and ECG holter	box of 1,200 (24 bags of 50)
33326	F2080P	Nonwoven	Ø 40 mm	Solid	Round	Snap	Pediatric	rest, long term, stress, ECG holter	box of 1,200 (24 bags of 50)
33327	F2080	Nonwoven	Ø 50 mm	Solid	Round	Snap	Adult	rest, long term, stress, ECG holter	box of 1,200 (24 bags of 50)
33341	F3001ECG	Gelled	34x23 mm	Solid	TAB	Alligator	Adult	rest ECG	bag of 100
33346	F3001ECG	Gelled	34x23 mm	Solid	TAB	Alligator	Adult	rest ECG	box of 5,000 (50 bags of 100)
33343	F3001NECG	Gelled	34x13 mm	Solid	TAB	Alligator	Pediatric, neonatal	rest ECG	box of 520 (13 bags of 40)
33372	F9047/4FL	PE-Foam	28x44 mm	Liquid	Rectang.	Plug 4 mm	Adult	rest ECG	box of 1,000 (20 bags of 50)
33373	F9053N	Foam	Ø 23x30 mm	Solid	Oval	4 mm socket	Pediatric, neonatal	monitoring ECG recording	box of 150 (5 bags of 30)
33377	F9058N	Foam	Ø 23x30 mm	Solid	Oval	1.5 mm socket	Pediatric, neonatal	long monitoring ECG recording	box of 150 (50 bags of 3)

*Radiolucent and amagnetic

ADAPTORS FOR ECG ELECTRODES

Male snap 4 mm

Alligator connector

Alligator connector

4 mm socket

Female 4 mm pin

Ø 4 mm socket

Front connection to TAB electrodes

4 mm socket

DOUBLE CONNECTION: SNAP OR TAB ELECTRODES

4 mm socket

Bottom connection to snap electrodes

LIMB CLAMP AND CHEST ELECTRODES

Silver coated plastic material Conform to EN 60601-2-25

Universal connection: hole with screw (4 mm) and clip

Universal connection - hole with screw (4 mm) and clip

Silver coated plastic material conform to EN 60601-2-25

30x130 mm adult clamp

20x90 mm pediatric clamp

PVC blue balls Ø 33 mm

PVC balls Ø 25 mm

Electrode Ø 14 mm

4 mm elastic connector

4 mm metal connector hole with screw and clips

GIMA code	UNIVERSAL REUSABLE ECG ELECTRODES LATEX FREE - NICKEL FREE	Minimum order
33363	Set of 4 limb clamp electrodes (4 colours) pediatric	1 set
33364	Set of 4 limb clamp electrodes (4 colours) adult	1 set
33368	Set of 6 chest electrodes Ø 15 mm - silver coated	1 set
33365	Set of 6 chest electrodes Ø 24 mm - silver coated	1 set
33369	Set of 6 chest electrodes Ø 30 mm - silver coated	1 set
33323	Set of 6 chest electr. Ultra Ø 24 mm - silver coated	1 set
33324	Set of 6 chest electr. Ultra Ø 30 mm - silver coated	1 set