

SURGICAL ABSORBABLE SUTURES - STERILE

MDR Certified

MONOGLYDE - monofilament absorbable suture poliglecaprone 25

Equivalent Ethicon brand: MONOCRYL

Technical details

- monofilament
- poliglecaprone 25
- synthetic absorbable suture
- sterilized by ethylene oxide
- shelf life of 3 years
- wound support: 21-28 days
- mass absorption: 90-120 days

Applications

- used for subcuticular repair and general soft-tissue approximation
- bowel anastomoses
- closure of the peritoneum
- closure of uterus and vaginal cuff
- not recommended for use in cardiovascular or neurological tissues, microsurgery or ophthalmic surgery
- since this is an absorbable suture, it should not be used for extended approximation of tissue (fascia)

Multilingual manual and box: GB, FR, IT, ES, PT, DE, GR.



MONOGLYDE - box of 12

GIMA code	Gauge	Metric	Colour	Circle	Needle profile	Needle length	Suture length
22901	0	3.5	Undyed	1/2	Round body	30 mm	70 cm
22902	2/0	3.0	Undyed	1/2	Round body	30 mm	70 cm
22903	3/0	2.0	Undyed	3/8	Slim blade cutting	25 mm	70 cm
22904	5/0	1.0	Undyed	3/8	Round body	19 mm	45 cm

TRUSYNTH FAST - braided & coated polyglactin 910 fast

Equivalent Ethicon brand: VICRYL RAPID

Technical details

- synthetic absorbable surgical suture
- braided and coated polyglactin 910
- sterilized by ethylene oxide
- shelf life of 5 years
- absorption period: 42 days

Applications

- general soft tissue approximation
- Multilingual manual and box: GB, FR, IT, ES, PT, DE, GR.



TRUSYNTH FAST - box of 12

GIMA code	Gauge	Metric	Colour	Circle	Needle profile	Needle length	Suture length
22910	2/0	3.0	Undyed	3/8	Reverse cutting	24 mm	70 cm
22911	3/0	2.0	Undyed	1/2	Round body	17 mm	45 cm
22912	3/0	2.0	Undyed	3/8	Reverse cutting	19 mm	45 cm
22913	3/0	2.0	Undyed	1/2	Round body	17 mm	70 cm
22914	3/0	2.0	Undyed	1/2	Round body	22 mm	70 cm
22915	3/0	2.0	Undyed	3/8	Reverse cutting	14 mm	75 cm
22916	3/0	2.0	Undyed	3/8	Reverse cutting	24 mm	70 cm
22917	4/0	1.5	Undyed	3/8	Reverse cutting	19 mm	45 cm
22918	4/0	1.5	Undyed	3/8	Reverse cutting	19 mm	70 cm
22919	5/0	1.0	Undyed	3/8	Reverse cutting	16 mm	70 cm
22920	6/0	0.7	Undyed	3/8	Reverse cutting	11 mm	45 cm

TRUGLYDE FAST - braided and coated polyglycolic acid fast undyed

Technical details

- braided and coated
- coated with calcium stearate and polycaprolactone
- sterilized by ethylene oxide
- shelf life of 5 years
- wound support: 7-10 days
- mass absorption: 40-42 days

Applications

- for episiotomy
- for circumcision
- closure of skin

Multilingual manual and box: GB, FR, IT, ES, PT, DE, GR.



TRUGLYDE FAST - box of 12

GIMA code	Gauge	Metric	Colour	Circle	Needle profile	Needle length	Suture length
22980	1	4.0	Undyed	1/2	Round body	40 mm	90 cm
22981	2/0	3.0	Undyed	1/2	Round body	35 mm	75 cm
22983	3/0	2.0	Undyed	3/8	Reverse cutting	24 mm	70 cm
22984	3/0	2.0	Undyed	3/8	Reverse cutting	20 mm	75 cm
22986	4/0	1.5	Undyed	3/8	Reverse cutting	16 mm	70 cm
22988	4/0	1.5	Undyed	3/8	Reverse cutting	20 mm	75 cm
22989	5/0	1.0	Undyed	1/2	Round body	17 mm	75 cm

PDSYNTH - monofilament absorbable suture polydioxanone

Equivalent Ethicon brand: PDS or PDS II

Technical details

- monofilament
- polydioxanone
- synthetic absorbable suture
- sterilized by ethylene oxide
- shelf life of 5 years
- wound support: 35-42 days
- mass absorption: 180-210 days

Applications

- for rectus sheath closure
- fascia closure and gastrointestinal tract
- pediatric cardiovascular surgery
- geriatric surgeries
- orthopedic surgery
- general surgery

Multilingual manual and box: GB, FR, IT, ES, PT, DE, GR.



PDSYNTH - box of 12

GIMA code	Gauge	Metric	Colour	Circle	Needle profile	Needle length	Suture length
22940	1	4.0	Violet	1/2	Round body	40 mm	90 cm
22941	0	3.5	Violet	1/2	Round body	30 mm	70 cm
22942	0	3.5	Violet	1/2	Round body	40 mm	70 cm

TRUSYNTH - braided & coated polyglactin 910

Equivalent Ethicon brand: VICRYL

Technical details

- braided and coated polyglactin 910
- coating of calcium stearate and polyglactin 910
- sterilized by ethylene oxide
- shelf life of 5 years
- wound support: 28-30 days
- mass absorption: 56-70 days

Applications

- ophthalmic procedures
- mucosal closures
- skin closure
- dental surgery
- orthopedic surgery
- general soft tissue approximation

Multilingual manual and box: GB, FR, IT, ES, PT, DE, GR.



TRUSYNTH - box of 12

GIMA code	Gauge	Metric	Colour	Circle	Needle profile	Needle length	Suture length
22950	1	4.0	Violet	1/2	Round body	40 mm	90 cm
22951	0	3.5	Violet	1/2	Round body	40 mm	90 cm
22952	2/0	3.0	Violet	1/2	Taper cutting	19 mm	75 cm
22953	2/0	3.0	Violet	1/2	Round body	30 mm	90 cm
22954	3/0	2.0	Violet	1/2	Taper cutting	19 mm	75 cm
22955	4/0	1.5	Violet	1/2	Taper cutting	19 mm	75 cm

TRUGLYDE - braided & coated polyglycolic acid

Technical details

- braided and coated polyglycolic acid
- coating of calcium stearate and polyglactin 310
- sterilized by ethylene oxide
- shelf life of 5 years
- wound support: 28-30 days
- mass absorption: 70-90 days

Applications

- for soft tissue approximation during general, gynecological and orthopedic surgeries

Multilingual manual and box: GB, FR, IT, ES, PT, DE, GR.



TRUGLYDE - box of 12

GIMA code	Gauge	Metric	Colour	Circle	Needle profile	Needle length	Suture length
22960	1	4.0	Violet	1/2	Round body	40 mm	90 cm
22961	0	3.5	Violet	1/2	Round body	40 mm	90 cm
22962	2/0	3.0	Violet	3/8	Reverse cutting	24 mm	70 cm
22963	2/0	3.0	Violet	1/2	Round body	30 mm	90 cm
22964	3/0	2.0	Violet	3/8	Reverse cutting	19 mm	45 cm
22965	3/0	2.0	Violet	3/8	Reverse cutting	19 mm	70 cm
22966	3/0	2.0	Undyed	3/8	Reverse cutting	19 mm	75 cm
22967	3/0	2.0	Undyed	1/2	Round body	22 mm	75 cm
22968	4/0	1.5	Violet	3/8	Reverse cutting	19 mm	45 cm
22969	4/0	1.5	Violet	3/8	Reverse cutting	19 mm	70 cm
22970	4/0	1.5	Undyed	3/8	Reverse cutting	16 mm	45 cm
22971	4/0	1.5	Undyed	1/2	Round body	17 mm	75 cm
22972	5/0	1.0	Undyed	1/2	Round body	17 mm	75 cm
22973	5/0	1.0	Undyed	3/8	Reverse cutting	13 mm	45 cm