

ESSEX FOLDING WHEELCHAIRS



Manual: GB, FR, IT, ES, PT, DE, GR, SE



- 43260 ESSEX WHEELCHAIR - 43 cm - black tissue
- 43261 ESSEX WHEELCHAIR - 46 cm - black tissue
- 43262 ESSEX WHEELCHAIR - 51 cm - black tissue

- blue powder coated steel frame
- rear solid PU wheels Ø 30.5 cm (12")
- front solid PU wheels Ø 20 cm (8")
- locking brakes on rear wheels
- front powder steel fork
- padded flip-up armrests
- detachable aluminium footrests with 6 height adjustable clips, flip-up plate, heel straps



ACCESSORIES FOR ESSEX

- 43263 Aluminium footplate - pair - spare
- 43264 Solid PU front wheel - spare
- 43265 Solid PU rear wheel - spare
- 43266 PVC arm pad - spare

Wheelchair	Seat			Wheels			Height		Depth		Weight	DIM		kg		
	W	D	H	Ø	Ø	H	W	D	W	D	(Kg)	W	D	NW	GW	
43260	58	25	43	30.5	20	50	91	96	40	41	110	76x25x93	16	18		
43261	61	25	46	30.5	20	50	91	96	40	41	110	76x25x93	16	18		
43262	66	25	51	30.5	20	50	91	96	40	41	110	76x25x93	16	18		

OXFORD FOLDING WHEELCHAIRS



Black nylon seat and back



- 43270 OXFORD WHEELCHAIR - 43 cm - black tissue
- 43271 OXFORD WHEELCHAIR - 46 cm - black tissue
- 43272 OXFORD WHEELCHAIR - 51 cm - black tissue

- blue powder coated steel frame
- pneumatic rear wheels Ø 61 cm (24") with chromed steel handrim and quick release axles
- front solid PU wheels Ø 20 cm (8")
- locking brakes on rear wheels
- front powder steel fork
- padded flip-up armrests
- detachable footrests, flip-up plate, with leg straps



Supplied with pump and manual: GB, FR, IT, ES, PT, DE, GR, SE.

- 43274 OXFORD PLUS WHEELCHAIR - 46 cm - black tissue
 - 43275 OXFORD PLUS WHEELCHAIR - 51 cm - black tissue
- Same features as above, but with solid PU rear wheels Ø 61 cm (24") and anti tipper system.

ACCESSORIES FOR OXFORD

- 43264 Solid PU front wheel - spare for 43270-5
- 43266 PVC arm pad - spare for 43270-5
- 43276 Plastic footplate - pair - spare
- 43277 Pneumatic rear wheel - spare for 43270-2
- 43278 Solid PU rear wheel - spare for 43274-5



Wheelchair	Seat			Wheels			Height		Depth		Weight	DIM		kg		
	W	D	H	Ø	Ø	H	W	D	W	D	(Kg)	W	D	NW	GW	
43270	60	26.5	43	61	20	50	91	103	40	44	120	82x27x93	18	20		
43271, 43274	63	26.5	46	61	20	50	91	103	40	44	120	82x27x93	18	20		
43272, 43275	68	26.5	51	61	20	50	91	103	40	44	120	82x27x93	18	20		

FOLDABLE COMMODOE WHEELCHAIR



43283 COMMODOE WHEELCHAIR

Chrome steel structure device ideal for patients with low mobility. It has fixed armrests and footrests, solid castors and sturdy rear wheels. U shape commode seat is extremely durable, suitable for indoor use. Supplied with toilet bowl.

Rear wheels diameter: 58 cm (22").

Front castors diameter: 20 cm (7").

- 27744 SOLID REAR WHEEL - spare
- 27766 PLASTIC WC BOWL Ø 26xh 15 cm - spare
- 43285 PLASTIC FOOT PLATE - pair - spare

Seat			Wheels			Height		Depth		Weight	DIM		kg	
W	D	H	Ø	Ø	H	W	D	W	D	(Kg)	W	D	NW	GW
64	44	58	20	50	87	105	37	100	92x29x89	100	19	21,3		

43283