



## TRAD STETHOSCOPES - GOOD VALUE FOR MONEY




**TRAD STETHOSCOPES: adult, pediatric**  
 Traditional aluminium stethoscopes, lightweight and reliable.  
 Chestpiece diaphragm Ø 43.5 mm (Ø 35 for pediatric model), bell Ø 30 mm (Ø 26.5 for pediatric model) external sizes.  
 Brass binaural, non chill ring.  
 Tubing: length 62 cm, int. Ø 4 mm, ext. Ø 8 mm.  
 Supplied in cardboard box. Manual and box in 11 languages: GB, FR, IT, ES, DE, PT, GR, PL, HU, SE, Arabic.



GIMA code	TRAD STETHOSCOPES - LATEX FREE	
32557	TRAD single head - pediatric - pink	
32559	TRAD single head - pediatric - light blue	
32555	TRAD double head - pediatric - black	
32560	TRAD single head - black	
32561	TRAD single head - red	
32562	TRAD single head - blue	
32563	TRAD single head - yellow	
32564	TRAD single head - purple	
32558	TRAD single head - dark green	
32565	TRAD double head - black	
32566	TRAD double head - red	
32567	TRAD double head - blue	
32568	TRAD double head - dark green	
32678	Diaphragm Ø 44 mm - spare	

## COLOURED TRAD STETHOSCOPES



**COLOURED TRAD STETHOSCOPES**  
 Coloured aluminium stethoscopes, lightweight and reliable.  
 Chestpiece diaphragm Ø 43.5 mm, bell Ø 30 mm external sizes.  
 Brass binaural.  
 Tubing: length 62 cm, int. Ø 4 mm, ext. Ø 8 mm.  
 Supplied in cardboard box.  
 Manual and box in 11 languages: GB, FR, IT, ES, DE, PT, GR, PL, HU, SE, Arabic.

GIMA code	COLOURED TRAD STETHOSCOPES - LATEX FREE	
51000	TRAD single head - black	
51001	TRAD single head - red	
51002	TRAD single head - blue	
51003	TRAD single head - yellow	
51004	TRAD single head - green	
51010	TRAD double head - black	
51011	TRAD double head - red	
51012	TRAD double head - blue	