

3M™ STERI-STRIP™ - REINFORCED ADHESIVE SKIN CLOSURE

Adhesive skin closures made of non-woven backing coated with a pressure-sensitive hypoallergenic adhesive. Sterile, breathable and comfortable to wear.



3M™ STERI-STRIP™ REINFORCED ADHESIVE SKIN CLOSURE
With polyester reinforcing filaments for higher tensile strength, reduced infection and scarring risk and better cosmetic outcomes.



34912



34912

STERI-STRIP™

GIMA code	3M code	3M™ REINFORCED STERI-STRIP™ - sterile - skin closure	Minimum order
34909	R1540	STERI-STRIP™ - 3x75 mm	50 bags of 5 strips
34910	R1541	STERI-STRIP™ - 6x75 mm	50 bags of 3 strips
34908	R1542	STERI-STRIP™ - 6x38 mm	50 bags of 6 strips
34911	R1546	STERI-STRIP™ - 6x100 mm	50 bags of 10 strips
34913	R1549	STERI-STRIP™ - 12x50 mm	50 bags of 6 strips
34912	R1547	STERI-STRIP™ - 12x100 mm	50 bags of 6 strips
34914	R1548	STERI-STRIP™ - 25x125 mm	25 bags of 4 strips

SUTURE STRIPS

• **34997 STERILE SUTURE STRIPS - 12 box of 8**
Sterile adhesive strips for cut and wound suture.
Cut to measure.
Dermatologically tested.
It contains 5 small strips (3x75 mm) and 3 medium strips (6.4x80 mm).
Multilingual box: GB, FR, IT, ES, DE.



34997

SOOTHING PATCH

• **34830 PREVIDOL SOOTHING PATCH**
9x14 cm - pouch of 5
Previdol is a Medical Device patch adjuvant in the treatment of pains in the joints and the muscular-skeletal system. It acts on the affected part, keeping the skin protected and helping to give relief in presence of pain.
Box: GB, IT.
To improve the cutaneous tolerability of the patch, some plant extracts like Arnica Montana, Devil's Claw and Boswellia were inserted, as well as some essential oils of Mint and Menthol to give a pleasant fragrance.



3M™ SKIN CARE SOLUTIONS - CAVILON™ ADVANCED

A complete range of barrier products to protect the skin from body fluids and friction. Cavilon™ products create on the skin a highly breathable, long lasting transparent coating to prevent maceration and trauma. Alcohol-free, hypoallergenic and non irritating, can be used on broken skin and under any tapes and dressing.

3M™ CAVILON™ CONTINENCE CARE WIPES
«All in 1» wipes for cleansing, moisturizing and protecting the skin.



34976

34974
3M™ CAVILON™ NO STING BARRIER FILM
Up to 72 hours of protection for peri-wound and peri-ostomy skin.



34975

3M™ CAVILON™ DURABLE BARRIER CREAM
Zinc-oxide free barrier cream for continence care. Apply twice a day.



GIMA code	3M code	3M™ CAVILON™	Minimum order
34974	3346E	Cavilon™ no sting barrier film - spray 28 ml	box of 12
34973	3346P	Cavilon™ no sting barrier film - spray 28 ml	box of 1
34977	3391G	Cavilon™ durable barrier cream 28 g	box of 1
34975	3392G	Cavilon™ durable barrier cream 92 g	1 dispenser
34976	9274	Cavilon™ continence care wipes	12 boxes of 8
34978	5051G4P	Cavilon™ high protection barrier film - liquid 0.7 ml	box of 4

3M™ CAVILON™ HIGH PROTECTION BARRIER FILM
Liquid solution with single-use applicator reduces the potential for cross-contamination and IAD pain creating a waterproof barrier film



34978

OMNISTRIP® - WOUND CLOSURE STRIPS

Omnistrip® is a wound closure strip for atraumatic closure of minor and surgical incisions in primary healing wounds. It protects the wound margins against tension and stabilizes fresh scars, once sutures or staples have been removed. Rounded edges, conformable and stretchable, with a high porosity, permeability and breathability and good adherence to the skin. Radiolucent, temperature resistant and durable. Made of a carded polyester-/polyamide skin-coloured non-woven, which is coated with an acrylic adhesive.
Class IIa, medical device. Multilingual manual and box: GB, FR, IT, ES, PT, DE, PL, NL, SE, DK, NO, FI, GR, CZ, SK, HU, RU, Arabic



GIMA code	Hartmann code	OMNISTRIP® WOUND CLOSURE STRIPS - STERILE	Minimum order
59300	540 681	OMNISTRIP® - 3x76 mm	box of 250 (50 bags of 5)
59301	540 682	OMNISTRIP® - 6x38 mm	box of 300 (50 bags of 6)
59302	540 683	OMNISTRIP® - 6x76 mm	box of 150 (50 bags of 3)
59303	540 684	OMNISTRIP® - 6x101 mm	box of 500 (50 bags of 2x5)
59304	540 685	OMNISTRIP® - 12x101 mm	box of 300 (50 bags of 2x3)
59305	540 686	OMNISTRIP® - 25x127 mm	box of 200 (50 bags of 2x2)



59300