

BED-TRAY - FOLDING LEGS

• **44751 BED TRAY - white**
 Plastic moulded bed-tray with containing edge.
 Chrome plated steel tube folding legs.
 Size: 58x39xh5 cm.
 Pack size: 65x45xh 8 cm.
 Weight: 1.9 kg



44751
 MADE IN ITALY

DINING TABLE

• **27449 DINING TABLE**
 Patient dining table with polished aluminium tubular structure. White laminated wooden top, with anti-shock rubber PVC band.
 Two rubber feet and two plastic wheels Ø 50 mm for easy mobility and comfortable use. Delivered in kit form.

Size: 110x42x93 cm

27449

MADE IN ITALY



OVERBED TABLES

Tray dimensions:
58x39xh 5 cm

MADE IN ITALY

Adjustable height:
75 to 114 cm

Base: 70 x 44 cm

27469

• **27469 OVERBED TABLE - moulded tray**
 Frame made of chrome plated steel tube. Serving tray made of thermoformed anti-shock plastic material with retaining edges and pre-formed bottle, glass and cutlery holder. Base with 4 Ø 50 mm castors. Carton size: 59x22x63 cm. Delivered in kit form.

Tabletop 60x40 cm

Adjustable height
61 to 94 cm

27475

4 castors,
2 with brake

MADE IN FRANCE

Tilting plate
0 to 90°

0-90°

0-15°
Tilting frame

• **27475 OVERBED DIFFUSION TABLE - white**
 Overbed table with high quality, laminated, fire resistant tabletop with an edging strip. Epoxy coated "C" shaped steel frame allows for closer access while sitting on a chair. Delivered ready to use and easy to fold flat (15 cm) for ease of storage. Weight: 8 kg. Max load: 10kg.

Adjustable height:
72 to 110 cm

Multilingual instructions:
GB, FR, IT, ES, PT,
DE, RO, SI.

27448

Set into the panel

Lock and adjustable panel

Adjustable height

• **27448 OVERBED TABLE**
 Overbed table with steel frame and high density board. Base with 4 Ø 50 mm (2.5") removable castors (2 with brake). Adjustable height and angle. Board size: 41x61 cm. Carton size: 65x7x67 cm. Delivered in kit form.

MASTER OVERBED TABLE

Height adjustment
(765-1050 mm)

Extra large ABS table top with cup socket

Air spring for height adjustment

27447

• **27447 MASTER OVERBED TABLE**
 High quality overbed table with adjustable height by gas spring system. Easy to clean table top made of ABS material, equipped with cup holder and high edges to protect objects from falling down. Anti-collision and durable round corner design. Bottom frame is sprayed with carbon steel. 4 revolving rubber castors Ø 50 mm, 2 with brake.
 Size: 890x440xh 765-1050 mm.

OLYMPUS BEDSIDE TABLES

43550

43551

43552

MADE IN ITALY

• **43550 OLYMPUS BEDSIDE TABLE**
 • **43551 OLYMPUS BEDSIDE TABLE with OVERBED TABLE**
 • **43552 OLYMPUS BEDSIDE TABLE with TILTING OVERBED TABLE**
 Double-sided bedside tables made of polyethylene. The inner compartment has a removable plastic shelf and a removable grid at the bottom to guarantee ventilation. Side panels include plastic towel holder, two cup and two bottle holders all equipped with retaining clamps. Drawer and doors made of moulded polyethylene plastic material with integrated handles. Height adjustable overbed tables (43551, 43552) with epoxy powder-coated steel frame with raised edges to avoid spills. 360° rotation and horizontal maneuverability of plate. Tilting top at different angles (43552 only). Four 55 mm Ø twin swivel castors.

TECHNICAL SPECIFICATIONS

Size: 61.5x48.5xh 75 cm (43550) 61.5x48.5xh 91 cm (43551-2) Overbed table (only 43551-2): size: 57x44 cm
 Weight: 14.5 kg (43550), 21.5 kg (43551-2) height: from 91 to 121 cm
 Drawer compartments size: small: 36x10xh 8.5 cm horizontal extension: 15 cm
 large: 36x18xh 8.5 cm
 Interior compartment size: 32.5x38xh 44 cm