


# GIMA DISPOSABLE STAINLESS STEEL OR CARBON STEEL BLADES


"GIMA" DISPOSABLE BLADES - sterile - dispenser box of 100  
Surgical blades are available in high quality carbon steel or stainless steel. Manufactured in accordance with requirements of British Standards BS 2982:1992, ISO 7740:1985, ISO 7153-1: 1991, EN 10088-3: 2005.  
Each blade is hallmarkd and individually packed in aluminium foil. Carbon Steel blades are VPI paper protected.  
Blade size is clearly identified on the reverse side of the foil peel-pack and on product unit boxes.  
Multilingual box and instructions: GB, FR, IT, ES, PT, DE, PL, HU, SK, SE, FI, DK, RO, BG, GR, HR, SI, Arabic. **MDR Certified.**



Carbon steel	Stainless steel	Diamantine stainless steel	GIMA STERILIZED BLADES	Minimum order
27000	27020	27010	Sterilized blades N. 10	box of 100
27001	27021	27011	Sterilized blades N. 11	box of 100
27002	27022	-	Sterilized blades N. 12	box of 100
27003	27023	27013	Sterilized blades N. 15	box of 100
27009	-	-	Sterilized blades N. 15c	box of 100
27004	27024	27014	Sterilized blades N. 20	box of 100
27005	27025	27015	Sterilized blades N. 21	box of 100
27006	27026	27016	Sterilized blades N. 22	box of 100
27007	27027	27017	Sterilized blades N. 23	box of 100
27008	27028	27018	Sterilized blades N. 24	box of 100



GIMA code	SCALPEL HANDLES - STAINLESS STEEL
26913	Scalpel Handle N° 3 for blades 10-15
26914	Scalpel Handle N° 4 for blades 20-25




**DIAMANTINE™ EDGE HIGH PERFORMANCE BLADES**  
High quality, super sharp blades with zirconium nitride coating fused on the cutting edge to satisfy the most demanding surgeons requirements. Durable cutting edge, that stays sharp incision after incision. The diamond like material creates a very hard edge that will not dull even under the most demanding cutting needs. Bio-compatible and highly corrosion resistant. Compatible with autoclavable scalpel handles N. 3 and N. 4. Stainless steel, sterile. **MDR Certified.**



# GIMA DISPOSABLE STAINLESS STEEL SCALPELS

"GIMA" DISPOSABLE SCALPELS - sterile - box of 10  
S/S blade fitted into an ABS plastic handle. Individually packed in sterilized envelope. Safety blade guard. **MDR Certified.**


"GIMA" SAFETY PROTECTIVE SHIELD SCALPELS - sterile  
Protective scalpels with a clear view shield that protects against injuries before use, between steps of procedure and after disposal. Latex and phthalates free. **MDR Certified**



Multilingual box:  
GB, FR, IT, ES, PT, DE

27040-27048

**PROTECTION GUARD**



Graduated ABS body



Perfectly balanced



Sonic-welded sandwiched blade




Anti skid blade guard

27060-27068


**STANDARD**

Protection guard	Standard	Safety	Diamantine high performance	GIMA STERILIZED DISPOSABLE STAINLESS STEEL SCALPELS (box of 10 pcs.)
27040	27060	27070	22460	Disposable scalpels - blade 10
27041	27061	27071	22461	Disposable scalpels - blade 11
27042	27062	-	-	Disposable scalpels - blade 12
27043	27063	27073	22463	Disposable scalpels - blade 15
27049	-	-	-	Disposable scalpels - blade 15c
27044	27064	27074	22464	Disposable scalpels - blade 20
27045	27065	27075	22465	Disposable scalpels - blade 21
27046	27066	27076	22466	Disposable scalpels - blade 22
27047	27067	27077	22467	Disposable scalpels - blade 23
27048	27068	27078	22468	Disposable scalpels - blade 24

• 27058 STITCH CUTTER - sterile - box of 10



27058 Disposable stainless steel stitch cutter



Clear view protective shield: blade clearly visible

Single handed action, easy to use for both right and left hand users

Complies with US OSHA standard 29 CFR 1910.1030 and European Directive 2010/32/EU

**SAFETY SCALPELS**


Stainless steel blade  
Stead fast disposable lock, audible, visual and tactile indication  
Blade size colour coded button  
Graduated handle (cm)

27070-27078

Multilingual box and instructions: GB, FR, IT, ES, PT, DE, NL, SI\*, SK\*, DK\*, LT\*, LV\*, EE\*, RU, RO, PL, BG, CZ, HU, SE, FI, GR, HR, Arabic. \* Only for standard scalpels

**DIAMANTINE HIGH PERFORMANCE SCALPELS - sterile**  
High-Performance Edge blade for extreme and durable sharpness  
- sensory grip handle enhances cutting action and control of the line of incision  
- each handle has a graduated 6 cm ruler to help determine incision length (indicative only)  
- protective sheath covers blade to help prevent accidental cuts  
- latex and phthalates free  
**MDR Certified**

**DIAMANTINE HIGH PERFORMANCE SCALPELS**



Stainless steel Diamantine Edge high performance coated blade. Hardness over 2,500 Vickers

22460-22468

Size is clearly marked on blade and handle for easy identification