

COUCH PAPER ROLLS

COUCH PAPER ROLLS


A complete range of hygienic and resistant couch paper manufactured in Italy in pure cellucotton or cellucotton/polythene. Each roll is individually packed in plastic bag for best protection. Atoxic. Available with standard 38 or 40 cm perforation or 190 cm perforation (full person size).

• 27418 COUCH ROLL HOLDER

Adjustable roll holder suitable for round (Ø 27-30 mm) or square (27-30 mm x 27-30 mm) metal couch leg.

MADE IN ITALY

OTHER SIZES OF COUCH ROLL PAPERS
AVAILABLE ON REQUEST



COUCH PAPER ROLLS

MATERIAL	Cellucotton embossed		Microembossed glued						Pure wadding of cellulose ply coated with a polythene protection					
	GIMA code	27410	27411	27412	27424	27427	27428	27432	27419	27414	27413	27415	27395	27396
Colour	white	white	white	white	white	white	white	white	white	light blue	white	white	light blue	white
Ply	1 ply	2 ply	2 ply	2 ply	2 ply	2 ply	2 ply	2 ply	2 ply	12 g/m² polyth.	12 g/m² polyth.	8 g/m² polyth.	12	7
Weight	30 g/m²	2x20 g/m²	2x17 g/m²	2x17 g/m²	2x17 g/m²	2x17 g/m²	2x17 g/m²	2x17 g/m²	2x17 g/m²	25/27 g/m²,	25/27 g/m²,	26/28 g/m²	23 g/m²	22 g/m²
Length	95 m	47.5 m	46 m	50 m	100 m	50 m	50 m	80 m	80 m	50 m	50 m	50 m	68 m	100 m
Height	50 cm	59 cm	50 cm	50 cm	50 cm	59 cm	59 cm	59 cm	59 cm	50 cm	50 cm	60 cm	50 cm	60 cm
N. of sheet	250	125	135	132	263	132	210	42	132	132	132	178	263	
Perforation	38 cm	38 cm	34 cm	38 cm	38 cm	38 cm	38 cm	190 cm	38 cm	38 cm	38 cm	38 cm	38 cm	38 cm
Rolls per box	6	6	9	9	6	9	6	6	6	6	6	6	6	6

SUITCASE COUCH

MADE IN ITALY



Couch size:
187.5x64xh 71 cm (open)
93x68xh 14 cm (folded)

• **27407 SUITCASE COUCH - black**
Chromed steel and folding legs.
Upholstered covering in washable material. Resin padded.
Max load: 100 kg.
Weight: 16.5 kg
Supplied with a velvet bag for easy transport.
Packing: 95x18x73.5 cm
Manual: GB, IT, NL.

ALUMINIUM SUITCASE COUCH

MADE IN ITALY



Manual: GB, FR, IT, ES, PT, DE, GR.



27406

• **27406 ALUMINIUM SUITCASE COUCH with hole**
Light, foldable, two sections, suitcase couch with aluminium frame. Central double legs for improved safety and block system on the hinges.
Head rest with mouth/nose hole.
Upholstered and covered with simulated leather.
Max load: 100 kg.
Weight: 13 kg.
Packing: 94x21.5x66.5 cm