


# NEWBORN TO OBESE CUFFS FOR ANEROID SPHYGMOMANOMETER



## CALIBRATED NYLON CUFFS

Calibrated nylon cuffs with TPU bladder and PVC tube.


CALIBRATED CUFFS - latex free - nylon					
1 tube*	2 tubes**	Colour	Model	Size cm	Arm circumference
-	32839	Black	Pediatric	35x11	19 to 26 cm
32846	32847	Blue	Pediatric	35x11	19 to 26 cm
32842	32843	Blue	Adult	49x15	26 to 40 cm
32844	-	Red	Adult	49x15	26 to 40 cm
32845	32850	Black	Adult	49x15	26 to 40 cm
32848	32849	Blue	Adult large	62x17.5	34 to 51 cm
-	32841	Black	Adult large	62x17.5	34 to 51 cm



## PREMIUM CUFFS

Made of high quality nylon 210D.  
Durability test up to 10,000 cycles.

PREMIUM CUFFS dark grey - nylon - TPU bladder - velcro fastening				
1 tube*	2 tubes**	Model	Size cm	Arm circumference
32817	32818	Newborn	5.5x18.5	9.5 to 13.5 cm
32886	32887	Infant	7.5x26	13 to 19 cm
32888	32889	Child	11x34.5	19 to 27 cm
32830	32832	Adult	54x15	25.4 to 40.6 cm
32834	32835	Adult large	61x17.5	34.3 to 50.8 cm
32837	32838	Thigh	82x21.5	40.6 to 66 cm
ADULT PREMIUM CUFFS coloured - nylon - TPU bladder - velcro fastening				
49750	-	Adult - royal blue	54x15	25.4 to 40.6 cm
49751	-	Adult - teal	54x15	25.4 to 40.6 cm
49752	-	Adult - ceil blue	54x15	25.4 to 40.6 cm
49753	-	Adult - pink	54x15	25.4 to 40.6 cm
ADULT CUFFS grey/green - cotton - TPU bladder - velcro fastening				
32829	-	Adult	54x15	25.4 to 40.6 cm




## BLADDERLESS NYLON CUFFS

- reusable above 10,000 times
- prevent cross-contaminations
- washable up to 60°C
- soft feeling on the skin
- disinfectable with alcohol/ethanol
- according to ISO 105 E01

BLADDERLESS ONE PIECE CUFF - latex free				
1 tube*	2 tubes**	Model	Size	Arm circumference
32910	32913	Children	10x34 cm	19 to 27 cm
32911	32914	Adult	14.5x54 cm	26 to 41 cm
32912	32915	Adult large	16x63 cm	34.3 to 51 cm
32919	32890	Thigh	20x83 cm	40.6 to 66.5 cm

\*One tube cuff for Sirio, London, Boston, Roma 1T, Tokyo, Madrid, Ri-san, Precisa, Minimus...  
\*\*Two tubes cuff for Minor, Yton line, Roma 2T, Dallas, Big Ben, Empire, Dayton...

# SPHYGMOMANOMETERS SPARE PARTS AND ACCESSORIES



GIMA code	SPHYGMOMANOMETERS SPARE PARTS	Minimum order
32865	Coiled tubing - black (42-46 spirals) - 39.5 cm	1
32808	Nylon velcro cuff (no bladder) 54 x15 cm*	1
49700	PVC Bladder - 1 tube - child	box of 5
49702	PVC Bladder - 1 tube - adult	box of 5
49703	PVC Bladder - 2 tubes - adult	box of 5
49706	TPU Bladder - 1 tube - adult	box of 5
49707	TPU Bladder - 2 tubes - adult	box of 5
32870	Male-Female metal connector	box 10
32871	Screw Type metal connector	box 10
32872	Bayonet Type metal connector	box 10
32873	Plastic connector	box 100

\*Can be used with 32874 or 32875

GIMA code	SPHYGMOMANOMETERS SPARE PARTS	Minimum order
32854	Yota PVC bulb - black - latex free	1
32851	Yota PVC bulb - green	1
32863	PVC bulb - black	box of 10
32855	Yota latex bulb - black	1
32860	Superior PVC bulb - green - latex free*	1
32859	Superior PVC bulb - black*	1
32852	Yota chrome-plated air release valve	1
32861	Superior chrome plated air valve	1
32853	Yota inflation system	1
32862	Superior inflation system - latex free*	1
32863	Deluxe inflation system	1

\*Made in Germany