

GIMACARE MULTIPARAMETER MONITOR - 7 PARAMETERS

• 24128 GIMACARE MULTI PARAMETER MONITOR

- 7 parameters - GB, FR, IT, ES

7 in 1 multiparameter diagnostic device for the periodic self-monitoring of various risk factors: blood glucose, β -ketone, lactate, total cholesterol, triglycerides, uric acid and hemoglobin. Ideal for both professional and homecare use. Results in 5 seconds.

Automatic strip identification system

Automatic strip recognition technology in one meter as guide to the right test

Bluetooth connected device

Data transmission thanks to Bluetooth connection with the development of miniaturization and wireless communication. Free connection to ProCheck App

3 coloured indication lights

1 colour indicator for correct strip insertion

- blue: strip guiding light

2 different colour indicators to understand whether results are low, high or in range

- green: within target range

- red: out of target range (low or high)

Additional features:

- automatic quality control

- before/after meal presettable

- strip ejector

- strip expiry date reminder

- auto sample loading detection

- auto electrode insertion detection

Supplied with lancing device, 10 sterile lancets, user guide, carrying case and user manual in GB, FR, IT, ES.

Strips not included.



7 PARAMETERS




GIMA code	STRIPS FOR GIMACARE - IN VIAL	Minimum order
24132	Glucose strips	vial of 50
24133	Ketone strips	vial of 25
24134	Lactate strips	vial of 25
24135	Cholesterol strips	vial of 25
24136	Uric acid strips	vial of 25
24137	Hemoglobin strips	vial of 25
24138	Triglycerides strips	vial of 25

Multilingual box and manual: GB, FR, IT, ES *Manual on request

GIMA code	CONTROL SOLUTIONS FOR GIMACARE	Minimum order
24139	Glucose control solution	1
24140	Ketone control solution	1
24141	Lactate control solution	1
24142	Cholesterol control solution	1
24143	Uric acid control solution	1
24144	Hemoglobin control solution	1

LANCETS		Minimum order
23916	Lancets 28G - sterile	bag of 100
23918	Lancets 30G - sterile	bag of 100
23921	Lancets 32G - sterile	bag of 100
32200	AAA alkaline battery	box of 4



TECHNICAL SPECIFICATIONS

Ketone warning: yes
Day average: 7, 14, 21, 28, 60, 90 days for blood glucose
Alarm: 4 daily alarms
Memory: 1000 measurements
Connection: Bluetooth
Power supply: 2 x 1.5 V AAA batteries
Size - weight: 102.5x56.9x21.8 mm - 57 g (without batteries)

PARAMETERS	MEASUREMENT RANGE	Sample size	Reaction time
Blood Glucose	10-800 mg/dL - 0.56-44.4 mmol/L	0.5 μ l	5 seconds
Ketone	0.1 - 8.0 mmol/L	0.8 μ l	10 seconds
Lactate	0.3 - 22 mmol/L	0.8 μ l	5 seconds
Cholesterol	100-400 mg/dL - 2.5-10.3 mmol/L	3.0 μ l	60 seconds
Triglycerides	70-600 mg/dL	3.0 μ l	60 seconds
Uric acid	3-20 mg/dL - 178-1190 μ mol/L	0.5 μ l	15 seconds
Hemoglobin	6.8 - 24 g/dL - 4.22-14.89 mmol/L	0.5 μ l	10 seconds

DRUG OF ABUSE, URINE TESTS FOR PROFESSIONAL USE

• 24553 MULTI DRUG PANEL TESTS - 6 drugs - urine

• 24554 MULTI DRUG PANEL TESTS - 8 drugs - urine

• 24556 MULTI DRUG PANEL TESTS - 10 drugs - urine

• 24557 MULTI DRUG PANEL TESTS - 12 drugs - urine

• 24559 MULTI DRUG PANEL TESTS - 16 drugs - urine

The Multi-Drug Rapid Test Panel is a rapid chromatographic immunoassay for the qualitative detection of multiple drugs and drug metabolites in human urine. Accurate results in 5 minutes. Multilingual manual and box: GB, FR, IT, ES, PT, DE, SE, PL, HU, RO, GR, Arabic.



SCREENING TEST PANEL					
Drug	Cut off ng/ml				
	24553	24554	24556	24557	24559
MOP	-	-	300	300	300
COC	150	150	300	150	300
THC	50	50	50	50	50
AMP	500	500	1000	500	1000
MDMA	-	500	-	500	500
MTD	-	-	300	300	300
MET	500	500	1000	500	1000
BZO	-	-	300	300	300
BAR	-	-	300	300	-
PCP	25	25	25	25	25
OPI	2000	2000	-	2000	-
TCA	-	-	1000	1000	-
6-AM	-	10	-	-	-
OXY	-	-	-	-	100
EDDP	-	-	-	-	300
BUP	-	-	-	-	-
TML	-	-	-	-	100
COT	-	-	-	-	200
K2	-	-	-	-	50
KET	-	-	-	-	1000
MDPV	-	-	-	-	1000