

GIMACARE MULTIPARAMETER MONITOR - 7 PARAMETERS

• 24128 GIMACARE MULTI PARAMETER MONITOR

- 7 parameters - GB, FR, IT, ES

7 in 1 multiparameter diagnostic device for the periodic self-monitoring of various risk factors: blood glucose, β -ketone, lactate, total cholesterol, triglycerides, uric acid and hemoglobin. Ideal for both professional and homecare use. Results in 5 seconds.

Automatic strip identification system

Automatic strip recognition technology in one meter as guide to the right test

Bluetooth connected device

Data transmission thanks to Bluetooth connection with the development of miniaturization and wireless communication. Free connection to ProCheck App

3 coloured indication lights

1 colour indicator for correct strip insertion

- blue: strip guiding light

2 different colour indicators to understand whether results are low, high or in range

- green: within target range

- red: out of target range (low or high)

Additional features:

- automatic quality control

- before/after meal presettable

- strip ejector

- strip expiry date reminder

- auto sample loading detection


- auto electrode insertion detection

Supplied with lancing device, 10 sterile lancets, user guide, carrying case and user manual in GB, FR, IT, ES.

Strips not included.



7 PARAMETERS




GIMA code	STRIPS FOR GIMACARE - IN VIAL	Minimum order
24132	Glucose strips	vial of 50
24133	Ketone strips	vial of 25
24134	Lactate strips	vial of 25
24135	Cholesterol strips	vial of 25
24136	Uric acid strips	vial of 25
24137	Hemoglobin strips	vial of 25
24138	Triglycerides strips	vial of 25

Multilingual box and manual: GB, FR, IT, ES *Manual on request

GIMA code	CONTROL SOLUTIONS FOR GIMACARE	Minimum order
24139	Glucose control solution	1
24140	Ketone control solution	1
24141	Lactate control solution	1
24142	Cholesterol control solution	1
24143	Uric acid control solution	1
24144	Hemoglobin control solution	1

LANCETS		
23916	Lancets 28G - sterile	bag of 100
23918	Lancets 30G - sterile	bag of 100
23921	Lancets 32G - sterile	bag of 100
32200	AAA alkaline battery	box of 4



TECHNICAL SPECIFICATIONS

Ketone warning: yes

Day average: 7, 14, 21, 28, 60, 90 days for blood glucose

Alarm: 4 daily alarms

Memory: 1000 measurements

Connection: Bluetooth

Power supply: 2 x 1.5 V AAA batteries

Size - weight: 102.5x56.9x21.8 mm - 57 g (without batteries)

PARAMETERS	MEASUREMENT RANGE	Sample size	Reaction time
Blood Glucose	10-800 mg/dL - 0.56-44.4 mmol/L	0.5 μ l	5 seconds
Ketone	0.1 - 8.0 mmol/L	0.8 μ l	10 seconds
Lactate	0.3 - 22 mmol/L	0.8 μ l	5 seconds
Cholesterol	100-400 mg/dL - 2.5-10.3 mmol/L	3.0 μ l	60 seconds
Triglycerides	70-600 mg/dL	3.0 μ l	60 seconds
Uric acid	3-20 mg/dL - 178-1190 μ mol/L	0.5 μ l	15 seconds
Hemoglobin	6.8 - 24 g/dL - 4.22-14.89 mmol/L	0.5 μ l	10 seconds

DRUG OF ABUSE, URINE TESTS FOR PROFESSIONAL USE

• 24553 MULTI DRUG PANEL TESTS - 6 drugs - urine

• 24554 MULTI DRUG PANEL TESTS - 8 drugs - urine

• 24556 MULTI DRUG PANEL TESTS - 10 drugs - urine

• 24557 MULTI DRUG PANEL TESTS - 12 drugs - urine

• 24559 MULTI DRUG PANEL TESTS - 16 drugs - urine

The Multi-Drug Rapid Test Panel is a rapid chromatographic immunoassay for the qualitative detection of multiple drugs and drug metabolites in human urine. Accurate results in 5 minutes. Multilingual manual and box: GB, FR, IT, ES, PT, DE, SE, PL, HU, RO, GR, Arabic.



SCREENING TEST PANEL

Drug	Cut off ng/ml					Detection time
	24553	24554	24556	24557	24559	
MOP	-	-	300	300	300	1-3 days
COC	150	150	300	150	300	1-7 days
THC	50	50	50	50	50	Use: single: 3 days; moderate: 5-7 days daily: 10-15 days; long term: >30 days
AMP	500	500	1000	500	1000	2-4+ days
MDMA	-	500	-	500	500	1-3 days
MTD	-	-	300	300	300	1-7 days
MET	500	500	1000	500	1000	1-4+ days
BZO	-	-	300	300	300	1-3 days, 4-6 weeks Chronic use
BAR	-	-	300	300	-	3-10 days
PCP	25	25	25	25	25	2-30 days
OPI	2000	2000	-	2000	-	1-3 days
TCA	-	-	1000	1000	-	4-21+ days
6-AM	-	10	-	-	-	6-8 hours
OXY	-	-	-	-	100	1-3 days
EDDP	-	-	-	-	300	4-5 days
BUP	-	-	-	-	-	3-5+ days
TML	-	-	-	-	100	1-3 days
COT	-	-	-	-	200	2-4 days
K2	-	-	-	-	50	1-3 days (single dose only)
KET	-	-	-	-	1000	2-4 days
MDPV	-	-	-	-	1000	2-3 days