

3M™ TRANSPORE™ WHITE TAPE



3M™ TRANSPORE™ WHITE SURGICAL TAPE

Gentle tape with strong adhesion and easy bi-directional tear for dressing applications.

Hypoallergenic and breathable.

Good initial and long-term adhesion (72 hours), easy to handle with gloves, water-resistant.


Suggested applications:

- securing dressings to damp or dry skin
- taping fragile, at-risk skin
- anchoring tubes and devices
- when tape width must be customized
- scissors not needed, patch resectable with hands



GIMA code	3M code	3M™ TRANSPORE™ WHITE TAPE	Minimum order
34950	1534-0	TRANSPORE™ 13 mm x 9.14 m	box of 24
34951	1534-1	TRANSPORE™ 25 mm x 9.14 m	box of 12
34952	1534-2	TRANSPORE™ 51 mm x 9.14 m	box of 6

3M™ MICROFOAM™



3M™ MICROFOAM™ SURGICAL TAPE


Highly conformable elastic foam tape that stretches for compression applications or securing dressings on challenging areas. Multi-directional stretch to accommodate swelling and promote comfort. Gentle, secure adhesion for anatomical areas with irregular contours. Water-resistant. Hypoallergenic and not made with natural rubber latex. Quickly conforms to joints and skin folds.

GIMA code	3M code	3M™ MICROFOAM™	Minimum order
59710	1528-1	MICROFOAM™ 25 mm x 5 m	box of 12
59711	1528-2	MICROFOAM™ 50 mm x 5 m	box of 6

SILK PLASTER ROLL

- 34787 SILK PLASTER ROLL
10 m x 2.5 cm - box of 12
- 34788 SILK PLASTER ROLL
10 m x 5 cm - box of 6

Silk surgical tape, white colour.



HYPOR POR PLASTER ROLLS



- 34806 SKINTONE FABRIC**
- 34804 NON WOVEN - 5 m**
- 34798 NON WOVEN - 9.14 m**

Each roll is supplied with a plastic roll cover except code 34798. Multilingual box: GB, FR, IT, ES, DE.

PLASTER ROLLS

Skintone fabric	Hypoallergenic Non woven	Size	Minimum order
34803	34804	5 m x 2.5 cm	box of 12 rolls
-	34798	9.14 m x 2.5 cm	box of 12 rolls
34806	34807	5 m x 5 cm	box of 6 rolls
34799	Shop Display: for showing up to 88 boxes (32 x 25 xh 63 mm)		

NON WOVEN PAPER TAPE ROLL



Suitable for dressing application for sensitive skin areas of the body

GIMA code	NON WOVEN PAPER TAPE ROLL HYPOALLERGENIC ADHESIVE	Minimum order
34789	Tape roll h 25 mm x 5 m	box of 12 pcs
34793	Tape roll h 25 mm x 9.14 m	box of 12 pcs
34794	Tape roll h 50 mm x 9.14 m	box of 6 pcs

ADHESIVE PLASTERS



GIMA code	HYPOALLERGENIC ADHESIVE PLASTERS	Minimum order
34809	Water resistant - 19x72 mm	2,000 (100 boxes of 20)
34814	Water resistant - 4 mixed shapes*	2,000 (100 boxes of 20)
34817	Sensitive skin - 5 mixed shapes**	3,200 (80 boxes of 40)

*19x72, 25x72, 9x38 mm, Ø 22mm

**19x72, 25x72, 9x38, 38x38 mm, Ø 22mm

PLASTER DISPENSER



- 34826 PLASTER DISPENSER - with 2 refills (one each type)
- 34827 REFILL 42 TEXTILE PLASTERS 19x72 mm - box of 6
- 34828 REFILL 36 TEXTILE PLASTERS 25x72 mm - box of 6

Dispenser for adhesive plasters ready for use. It allows operator to extract plaster with one hand and apply directly on the wound. Transparent cover protects plaster and shows residual quantities. Lockable and with fixing wall kit included.