

CU Medical Systems, Inc. MADE IN KOREA

DEFIBRILLATOR
MODEL : CU-SP2

POWER :

Disposable Battery

≡ 12V, 4.2Ah

Rechargeable Battery

≡ 11.1V, 1.9Ah

IP55



2460



CU Medical Systems, Inc.

130-1, Donghwagongdan-ro, Munmak-eup,
Wonju-si, Gangwon-do, Korea

Tel: +82 33 747 7657 Fax: +82 33 747 7659

<http://www.cu911.com>

EC REP

Medical Device Safety Service

Schiffgraben 41, 30175 Hannover, Germany

Tel: +49 511 6262 8630 Fax: +49 511 6262 8633

SN



WARNING-ELECTRIC SHOCK HAZARD.

Do not remove cover(or back). Refer servicing to qualified service personnel.



WARNING-Possible explosion hazard if used in the presence of concentrated oxygen.