



OPEN HERE



STERILE URINARY CATHETER FOR SINGLE USE (Foley catheter) CATETERE URINARIO STERILE MONOUSO

WARNING: Do not use ointments or lubricants having a petroleum base. They will damage the catheter and cause the balloon to burst.

ATTENZIONE: Non utilizzare unguenti o lubrificanti a base di petrolio. Danneggeranno il catetere e faranno scoppiare il palloncino.

DIRECTIONS: After aseptic catheterization, inflate the retention balloon by inserting the syringe tip into the valve and introduce the necessary amount of sterile water.

INDICAZIONI: Dopo il cateterismo asettico gonfiare il palloncino di ritenzione inserendo la punta della siringa nella valvola e introdurre la quantità necessaria di acqua sterile.



Jiangsu Suyun Medical Materials Co., Ltd.
No.18 Jin Qiao Road Dapu Industrial Park 222002 Lianyungang, Jiangsu Province
PEOPLE'S REPUBLIC OF CHINA
Tel.: +86-518-85608151 Fax: +86-518-85466033 Made in China.
Shanghai International Holding Corp. GmbH (Europe)
Eiffestraße 80, 20537 Hamburg, GERMANY
Gima S.p.A.
Via Marconi, 1
20060 Gessate (MI) - Italy
www.gimaitaly.com

REV.2023-05

Type: 2-WAY-STANDARD TIP
Size: 4.0mm(12Fr/CH). 5-10cc/ml
REF S5122.12 (GIMA 22054)

LOT 000009
230909
280908



Type: 2-WAY-STANDARD TIP
Size: 4.0mm(12Fr/CH). 5-10cc/ml
REF S5122.12 (GIMA 22054)

LOT 000009
230909
280908



Type: 2-WAY-STANDARD TIP
Size: 5.3mm(16Fr/CH). 5-10cc/ml
REF S5122.16 (GIMA 22056)

LOT 000011
230909
280908



Type: 2-WAY-STANDARD TIP
Size: 6.0mm(18Fr/CH). 5-10cc/ml
REF S5122.18 (GIMA 22057)

LOT 000010
230909
280908



505.0 mm

70.0 mm

Type: 2-WAY-STANDARD TIP
Size: 4.0mm(12Fr/CH). 5-10cc/ml

REF S5122.12 (GIMA 22054)
LOT 000009

230909
280908

UDI (01)16933298298662
(10)000009
(11)230909
(17)280908

8 023279 220544

Type: 2-WAY-STANDARD TIP
Size: 5.3mm(16Fr/CH). 5-10cc/ml

REF S5122.16 (GIMA 22056)
LOT 000011

230909
280908

UDI (01)16933298298679
(10)000011
(11)230909
(17)280908

8 023279 220568

Type: 2-WAY-STANDARD TIP
Size: 6.0mm(18Fr/CH). 5-10cc/ml

REF S5122.18 (GIMA 22057)
LOT 000010

230909
280908

UDI (01)16933298298686
(10)000010
(11)230909
(17)280908

8 023279 220575

Rev.2022-05

Shanghai International Holding Corp. GmbH (Europe)
Eiffestraße 80, 20537 Hamburg, GERMANY

EC REP

STERILE URINARY CATHETER
FOR SINGLE USE
(Foley catheter)
CATETERE URINARIO STERILE MONOUSO

Gima S.p.A.
Via Marconi, 1
20060 Gessate (MI) – Italy
www.gimaitaly.com

Jiangsu Suyun Medical Materials Co., Ltd.
No.18 Jin Qiao Road Dapu Industrial Park 222002 Lianyungang, Jiangsu Province
PEOPLE'S REPUBLIC OF CHINA
Tel.: +86-518-85608151 Fax: +86-518-85466033
Made in China.

Medsuyun

STERILE URINARY CATHETER
FOR SINGLE USE
(Foley catheter)
CATETERE URINARIO STERILE MONOUSO

WARNING: Do not use ointments or lubricants having a petroleum base. They will damage the catheter and cause the balloon to burst.

ATTENZIONE: Non utilizzare unguenti o lubrificanti a base di petrolio. Danneggeranno il catetere e faranno scoppiare il palloncino.

DIRECTIONS: After aseptic catheterization, inflate the retention balloon by inserting the syringe tip into the valve and introduce the necessary amount of sterile water.

INDICAZIONI: Dopo il cateterismo asettico gonfiare il palloncino di ritenzione inserendo la punta della siringa nella valvola e introdurre la quantità necessaria di acqua sterile.

10PCS/BOX

STERILE EO MD 230909 280908

UDI (01)16933298298662
(10)000009
(11)230909
(17)280908

8 023279 220544

Medsuyun

STERILE URINARY CATHETER
FOR SINGLE USE
(Foley catheter)
CATETERE URINARIO STERILE MONOUSO

Gima S.p.A.
Via Marconi, 1
20060 Gessate (MI) – Italy
www.gimaitaly.com

Jiangsu Suyun Medical Materials Co., Ltd.
No.18 Jin Qiao Road Dapu Industrial Park 222002 Lianyungang, Jiangsu Province
PEOPLE'S REPUBLIC OF CHINA
Tel.: +86-518-85608151 Fax: +86-518-85466033 Made in China.

EC REP Shanghai International Holding Corp. GmbH (Europe)
Eiffestraße 80, 20537 Hamburg, GERMANY

Rev.2022-05