

**DEVICE LABEL**




ID label on the battery cover

(01)08052990321766  
(21)A23-B000001

UDI label (HRI format) on the back of the battery compartment  
VARIABLE with Serial Number

**PACKAGING LABEL**



write on the right side:

**MIR S.p.A.**  
Viale Luigi Schiavonetti  
00173 Rome (Italy)  
[www.spirometry.com](http://www.spirometry.com)



CE 0476
















write on the right side:

FCC ID: 2AAQS-ISP1507



consult instructions for use

|   |   |
|---|---|
|    | <p>applied part type</p>                                    |
| <p><b>IP 22</b></p>   | <p>degree of resistance to water and particles</p>          |
|    | <p>medical device symbol</p>                                |
|   | <p>do not expose to direct light</p>                        |
|  | <p>do not expose to atmospheric agents</p>                  |
|  | <p>Temperature range</p> <p>MIN: -25 °C<br/>MAX: 70 °C</p>  |
|  | <p>relative humidity range</p> <p>MIN: 10%<br/>MAX: 90%</p> |

|   |  |  |
|---|--|--|
|    | <p>sorted garbage</p>  |  |
| <p><b>Battery CR2032</b></p>  | <p>internal battery specification</p>  |  |
|    | <p>catalog reference<br/>write on the right side:<br/>91200000</p>                   |  |
|  <p>symbol indicating that the information next to it represents the UDI</p> |   | <p>Area of at least 3x4 cm where to place the 3 elements below</p>   |
|   |  | <p>put into production in the free rectangular area above</p> <p>UDI_DI (01) - FIXED</p> <p>SERIALE (21) - VARIABLE</p>  |
|   |   | <p>year of manufacture</p> <p>1. write on the right side: YYYY/MM</p> <p>2. place in the free rectangular area above</p> |
|   |  | <p>put into production in the free rectangular area above</p> <p>serial type:<br/><b>A23-ABNNNNN</b></p>                 |

ETICHETTA 01 PRINCIPALE DA APPORRE SOTTO A PACK - 174 x 55 mm

AREA UDI - 47,8 x 44 mm

MATRICE UDI - 42 x 13,65 mm

|             |                     |
|-------------|---------------------|
| matrice UDI | (01) 08052990321766 |
|             | (21) A23-B000001    |

55 mm

MIR S.p.A.  
viale Luigi Schiavonetti 270  
00173 Rome, Italy  
www.spirometry.com

FCC ID: 2AAQS-ISP1507  
Battery CR2032  
REF 91200000

UDI

IP 22 MD

10% 95% -25 °C 70 °C

BAR CODE - 42 x 10 mm



ANNO DI FABBRICAZIONE  
22 x 10 mm



174 mm