

LOT
CODE
UDI



8 023279 241341

GimaCare

- Usare solo con sistema di monitoraggio multifunzionale GIMACARE.
Strisce test, chip codice e istruzioni all'interno
- Use only with GIMACARE multi-functional monitoring system.
Test strips, code chip, and Instruction inside
- Utiliser uniquement avec système de surveillance multifonctionnel GIMACARE.
Bandelettes de test, puce de code et instructions à l'intérieur
- Usar únicamente con sistema de monitorización multifuncional GIMACARE.
Las tiras de prueba, el chip del código y las instrucciones están en el interior


Chemical components in sensor:
1. Lactate oxidase (Microorganism)
2. Mediator
3. Enzyme protector
4. Non-reactive ingredients

≥ 0.5 U
52%
6%
32%

www.gimaitaly.com



Campione di sangue
di piccole dimensioni 0,8 µ L
0.8 µ L Tiny Blood Sample
Échantillon de sang très petit 0,8 µ L
Muestra de sangre pequeña 0,8 µ L



Applicare una goccia di sangue qui
Apply a drop of blood here
Appliquer une goutte de sang ici
Aplique una gota de sangre aquí

GimaCare

LAC **Lactate
Test Strip**

Strisce Test
Lattato
Bandelettes
de Test Lactate
Tiras de Prueba
Lactato

25 PZ PCS
PCS PIEZAS

www.gimaitaly.com



GIMA 24134

TaiDoc Technology Corporation
B1-7F, No. 127, Wugong 2nd Rd., Wugu Dist.,
24888 New Taipei City, Taiwan
www.taidoc.com - Made in Taiwan

REF TD-4663

EC REP

CE
0123

MedNet EC-REP GmbH
Borkstraße 10, 48163 Münster, Germania
Importato da / Imported by / Importé par / Importado por:
Gima S.p.A.
Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com



Per automonitoraggio
For self-testing
Pour autosurveillance
Para automonitorización



Nota:
1. Non applicare sangue sulla parte superiore della striscia.
2. Qualsiasi area della striscia può essere toccata **senza alterare la lettura.**



Nota:
1. Do not apply blood on the top of the strip.
2. Any area of this strip can be touched and the **reading will not be affected.**

Nota:
1. Ne pas appliquer de sang sur le dessus de la bande.
2. Toute zone de la bandelette peut être touchée **sans altérer la lecture**

Nota:
1. No aplique sangre sobre la parte superior de la tira.
2. Cualquier área de la tira puede ser tocada **sin alterar la lectura**

509685 Rev. 0 - 10/24
Size 4x4,3x5,7 cm