

505.0 mm

70.0 mm



Medsuyun

STERILE URINARY CATHETER FOR SINGLE USE (Foley catheter) CATETERE URINARIO STERILE MONOUSO

WARNING: Do not use ointments or lubricants having a petroleum base. They will damage the catheter and cause the balloon to burst.

ATTENZIONE: Non utilizzare unguenti o lubrificanti a base di petrolio. Danneggeranno il catetere e faranno scoppiare il palloncino.

DIRECTIONS: After aseptic catheterization, inflate the retention balloon by inserting the syringe tip into the valve and introduce the necessary amount of sterile water.

INDICAZIONI: Dopo il cateterismo asettico gonfiare il palloncino di ritenzione inserendo la punta della siringa nella valvola e introdurre la quantità necessaria di acqua sterile.



Jiangsu Suyun Medical Materials Co., Ltd.
No.18 Jin Qiao Road Dapu Industrial Park 222002 Lianyungang, Jiangsu Province
PEOPLE'S REPUBLIC OF CHINA
Tel.: +86-518-85606151 Fax: +86-518-85466033 Made in China.

Shanghai International Holding Corp. GmbH (Europe)
Eiffelstraße 80, 20537 Hamburg, GERMANY

Gima S.p.A.
Via Marconi, 1
20060 Gessate (MI) - Italy
www.gimaitaly.com

REV.2023-05

Type: 2-WAY-STANDARD TIP
Size: 5.3mm(16Fr/CH). 30cc/ml

REF S5122.16 (GIMA 22066)

LOT 000013

230909

280908



(01)06933298298696
(10)000013
(11)230909
(17)280908

Type: 2-WAY-STANDARD TIP
Size: 5.3mm(16Fr/CH). 30cc/ml

REF S5122.16 (GIMA 22066)

LOT 000013

230909

280908

UDI



(01)06933298298696
(10)000013
(11)230909
(17)280908

Type: 2-WAY-STANDARD TIP
Size: 6.0mm(18Fr/CH). 30cc/ml

REF S5122.18 (GIMA 22067)

LOT 000015

230909

280908

UDI



(01)06933298298702
(10)000015
(11)230909
(17)280908

Type: 2-WAY-STANDARD TIP
Size: 8.0mm(24Fr/CH). 30cc/ml

REF S5122.24 (GIMA 22070)

LOT 000003

230909

280908



(01)06933298298719
(10)000003
(11)230909
(17)280908

Type: 2-WAY-STANDARD TIP
Size: 5.3mm(16Fr/CH), 30cc/ml

REF S5122.16 (GIMA 22066)
LOT 000013

230909
280908

UDI (01)16933298298716
(10)000013
(11)230909
(17)280908

8 023279 220667

Type: 2-WAY-STANDARD TIP
Size: 6.0mm(18Fr/CH), 30cc/ml

REF S5122.18 (GIMA 22067)
LOT 000015

230909
280908

UDI (01)16933298298709
(10)000015
(11)230909
(17)280908

8 023279 220674

Type: 2-WAY-STANDARD TIP
Size: 8.0mm(24Fr/CH), 30cc/ml

REF S5122.24 (GIMA 22070)
LOT 000003

230909
280908

UDI (01)16933298298716
(10)000003
(11)230909
(17)280908

8 023279 220704

EC REF
Shanghai International Holding Corp. GmbH (Europe)
Eiffestraße 80, 20537 Hamburg, GERMANY
Rev.2022-05

No.18 Jin Qiao Road Dapu Industrial Park 222002 Lianyungang,Jiangsu Province
PEOPLE'S REPUBLIC OF CHINA
Tel.:+86-518-85608151 Fax:+86-518-85466033
Made in China.

Jiangsu Suyun Medical Materials Co., Ltd

Gima S.p.A.
Via Marconi, 1
20060 Gessate (MI) – Italy
www.gimaitaly.com

STERILE URINARY CATHETER
(Foley catheter)
FOR SINGLE USE
CATETERE URINARIO STERILE MONOUSO



STERILE URINARY CATHETER
FOR SINGLE USE
(Foley catheter)
CATETERE URINARIO STERILE MONOUSO

WARNING: Do not use ointments or lubricants having a petroleum base. They will damage the catheter and cause the balloon to burst.

ATTENZIONE: Non utilizzare unguenti o lubrificanti a base di petrolio. Danneggeranno il catetere e faranno scoppiare il palloncino.

DIRECTIONS: After aseptic catheterization, inflate the retention balloon by inserting the syringe tip into the valve and introduce the necessary amount of sterile water.

INDICAZIONI: Dopo il cateterismo asettico gonfiare il palloncino di ritenzione inserendo la punta della siringa nella valvola e introdurre la quantità necessaria di acqua sterile.

STERILE EO MD 2 35°C +5°C 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19 20 21 22 23 24 25 26 27 28 29 30 31 32 33 34 35 36 37 38 39 40 41 42 43 44 45 46 47 48 49 50 51 52 53 54 55 56 57 58 59 60 61 62 63 64 65 66 67 68 69 70 71 72 73 74 75 76 77 78 79 80 81 82 83 84 85 86 87 88 89 90 91 92 93 94 95 96 97 98 99 100

50PCS/BOX

Type: 2-WAY-STANDARD TIP
Size: 5.3mm(16Fr/CH), 30cc/ml

REF S5122.16 (GIMA 22066)
LOT 000013

230909
280908

UDI (01)16933298298693
(10)000013
(11)230909
(17)280908

8 023279 220667



STERILE URINARY CATHETER
FOR SINGLE USE
(Foley catheter)
CATETERE URINARIO STERILE MONOUSO

Gima S.p.A.
Via Marconi, 1
20060 Gessate (MI) – Italy
www.gimaitaly.com

Jiangsu Suyun Medical Materials Co., Ltd.
No.18 Jin Qiao Road Dapu Industrial Park 222002 Lianyungang,Jiangsu Province
PEOPLE'S REPUBLIC OF CHINA
Tel.:+86-518-85608151 Fax:+86-518-85466033
Made in China.

EC REF
Shanghai International Holding Corp. GmbH (Europe)
Eiffestraße 80, 20537 Hamburg, GERMANY

Rev.2022-05