

工厂打印标签示意图

MD Full Face Mask/ Maschera facciale completa
Model/Modello: YF-01

REF MF01002

LOT 1231008 2023-10-17 2028-10-16

Suzhou Yuyue Medical Technology Co., Ltd.
No.9 Jinfeng Road, Suzhou Science & Technology Town,
215163 Suzhou, Jiangsu, PRC

EC REP Metrax GmbH
Rheinwaldstr. 22, 78628 Rottweil, Germany

UDI (01)06922166510589
(11)231017
(17)281016
(10)1231008

90% 10% 0123

MD Full Face Mask/ Maschera facciale completa
Model/Modello: YF-01

REF MF01003

LOT 1231008 2023-10-17 2028-10-16

Suzhou Yuyue Medical Technology Co., Ltd.
No.9 Jinfeng Road, Suzhou Science & Technology Town,
215163 Suzhou, Jiangsu, PRC

EC REP Metrax GmbH
Rheinwaldstr. 22, 78628 Rottweil, Germany

UDI (01)06922166510596
(11)231017
(17)281016
(10)1231008

90% 10% 0123

技术要求:

- 印刷内容以电子矢量图档为准;
- 标签打印清晰、无缺印、重印,切边整齐完整;
- 产品名称:根据实际生产填写,

口鼻罩(YF/RT系列)填写: Full Face Mask,
鼻罩(YN系列)填写: Nasal Mask
鼻枕面罩(YP系列)填写: Nasal pillows Mask

- 型号:根据实际生产打印产品型号,注明示例:YP-01;
- 规格:根据实际需求打印产品规格
L\M\S(规格外的△和L\M\S一起打印),组合销售的YP
不打印规格。

- 处为生产批次,注明示例: 1231008
生产批次为2023年10月第八批次;

- 处为生产日期,根据实际生产填写。

- 处为使用期限,使用期限为5年,注明示例:

注明示例:2023-10-15,使用期限为2023年10月15日(对应的生产日期为2023年10月16日)

9. UDI和REF号详见面罩REF清单

10. 标签跳距2.5mm,左右留边各2.5mm,公差为2.5mm+0.5mm/-0mm。

标记	处数	更改文件号	签 字	日期	
设 计			视 角		
制 图			工 艺		
审 核			版 本	A	
批 准			日 期		

面罩产品标签-英意双语-A

PET+可打印哑膜

YY05-C-B28

图样标记 数量 重量 比例

S 1 1:1

共 2 张 第 1 张

yuwell 鱼跃
让科技律动生命
苏州鱼跃医疗科技有限公司