



GIMA

Via Marconi, 1 - 20060 Gessate (MI) Italia

Web Site: www.gimaitaly.com



Rev.2-09.22

**LETTINO DA VISITA ELETTRICO AD ALTEZZA
REGOLABILE con barra perimetrale - blu**
**ELECTRIC HEIGHT ADJUSTABLE TREATMENT
TABLE with footbar - blue**

REF 44520

LOT XXXXXXXXX

SN XXXX



aaaa/mm



Gima S.p.A.

Via Marconi, 1 - 20060 Gessate (MI) Italy

gima@gimaitaly.com - export@gimaitaly.com

www.gimaitaly.com

Made in Italy

100-240V 1.5A 50/60Hz

2 min ON / 18min OFF (10%)

IP44 MD

