



OPEN HERE



STERILE URINARY CATHETER FOR SINGLE USE (Foley catheter)

CATETERE URINARIO STERILE MONOUSO

WARNING: Do not use ointments or lubricants having a petroleum base. They will damage the catheter and cause the balloon to burst.

ATTENZIONE: Non utilizzare unguenti o lubrificanti a base di petrolio. Danneggeranno il catetere e faranno scoppiare il palloncino.

DIRECTIONS: After aseptic catheterization, inflate the retention balloon by inserting the syringe tip into the valve and introduce the necessary amount of sterile water.

INDICAZIONI: Dopo il cateterismo asettico gonfiare il palloncino di ritenzione inserendo la punta della siringa nella valvola e introdurre la quantità necessaria di acqua sterile.



Jiangsu Suyun Medical Materials Co., Ltd.
No.18 Jin Qiao Road Dapu Industrial Park 222002 Lianyungang, Jiangsu Province
PEOPLE'S REPUBLIC OF CHINA
Tel.: +86-518-85608151 Fax: +86-518-85466033 Made in China.
Shanghai International Holding Corp. GmbH (Europe)
Eiffestraße 80, 20537 Hamburg, GERMANY

Gima S.p.A.
Via Marconi, 1
20060 Gessate (MI) - Italy
www.gimaitaly.com

REV.2023-05

Type: 3-WAY-STANDARD TIP
Size: 6.7mm(20Fr/CH). 30cc/ml

REF S5123.20 (GIMA 22078)

LOT 000009

230909

280908



(01)06933298298726
(10)000009
(11)230909
(17)280908

505.0 mm

70.0 mm

Type: 2-WAY-STANDARD TIP
Size: 5.3mm(16Fr/CH), 30cc/ml
REF S5122.16 (GIMA 22066)
LOT 000013
230909
280908

UDI
(01)16933298298716
(10)000013
(11)230909
(17)280908

8 023279 220667

Type: 2-WAY-STANDARD TIP
Size: 6.0mm(18Fr/CH), 30cc/ml
REF S5122.18 (GIMA 22067)
LOT 000015
230909
280908

UDI
(01)16933298298709
(10)000015
(11)230909
(17)280908

8 023279 220674

Type: 2-WAY-STANDARD TIP
Size: 8.0mm(24Fr/CH), 30cc/ml
REF S5122.24 (GIMA 22070)
LOT 000003
230909
280908

UDI
(01)16933298298716
(10)000003
(11)230909
(17)280908

8 023279 220704

EC REF
Shanghai International Holding Corp. GmbH (Europe)
Eiffestraße 80, 20537 Hamburg, GERMANY
Rev.2022-05

Jiangsu Suyun Medical Materials Co., Ltd.
No.18 Jin Qiao Road Dapu Industrial Park 222002 Lianyungang, Jiangsu Province
PEOPLE'S REPUBLIC OF CHINA
Tel.: +86-518-85608151 Fax: +86-518-85466033
Made in China.

Gima S.p.A.
Via Marconi, 1
20060 Gessate (MI) – Italy
www.gimaitaly.com

STERILE URINARY CATHETER
(Foley catheter)
FOR SINGLE USE
CATETERE URINARIO STERILE MONOUSO



STERILE URINARY CATHETER
FOR SINGLE USE
(Foley catheter)
CATETERE URINARIO STERILE MONOUSO

WARNING: Do not use ointments or lubricants having a petroleum base. They will damage the catheter and cause the balloon to burst.

ATTENZIONE: Non utilizzare unguenti o lubrificanti a base di petrolio. Danneggeranno il catetere e faranno scoppiare il palloncino.

DIRECTIONS: After aseptic catheterization, inflate the retention balloon by inserting the syringe tip into the valve and introduce the necessary amount of sterile water.

INDICAZIONI: Dopo il cateterismo asettico gonfiare il palloncino di ritenzione inserendo la punta della siringa nella valvola e introdurre la quantità necessaria di acqua sterile.



50PCS/BOX

Type: 2-WAY-STANDARD TIP
Size: 5.3mm(16Fr/CH), 30cc/ml
REF S5122.16 (GIMA 22066)
LOT 000013
230909
280908

UDI
(01)16933298298693
(10)000013
(11)230909
(17)280908

8 023279 220667



STERILE URINARY CATHETER
FOR SINGLE USE
(Foley catheter)
CATETERE URINARIO STERILE MONOUSO

Gima S.p.A.
Via Marconi, 1
20060 Gessate (MI) – Italy
www.gimaitaly.com

Jiangsu Suyun Medical Materials Co., Ltd.
No.18 Jin Qiao Road Dapu Industrial Park 222002 Lianyungang, Jiangsu Province
PEOPLE'S REPUBLIC OF CHINA
Tel.: +86-518-85608151 Fax: +86-518-85466033
Made in China.

EC REF
Shanghai International Holding Corp. GmbH (Europe)
Eiffestraße 80, 20537 Hamburg, GERMANY

Rev.2022-05