

EXCLUSIVELY MANUFACTURED FOR



GIMA

GIMA Spa - Via Marconi, 1
20060 Gessate (MI) - Italy

CUT

FREQUENCY

160W → 250 Ω 600 KHz

BLEND

FREQUENCY

120W → 200 Ω 600 KHz

FORCED COAG

FREQUENCY

100W → 150 Ω 600 KHz

SOFT COAG

FREQUENCY

80W → 100 Ω 600 KHz

BIPOLAR

FREQUENCY

60W → 100 Ω 600 KHz

SN

394



LED SpA
Via Selciatella 40
APRILIA (LT) - ITALY -

CLASS

MODEL

I CF

DIATERMO
MB 160 D

DUTY CYCLE

ON: 10 s / OFF: 30 s

WARNING
NO OPERATOR SERVICEABLE
PARTS INSIDE, REFER SERVICING
TO QUALIFIED PERSONNEL

ENSURE THAT THE CORRECT
VOLTAGE IS SELECTED BEFORE
CONNECTING TO THE MAINS
POWER SUPPLY



115 V~ - 50/60 Hz
350 VA - FUSE 2 X T6.3AL 250V

230 V~ - 50/60 Hz
350 VA - FUSE 2 X T3.15AL 250V

MADE IN
ITALY