

WARNING

NO OPERATOR
SERVICEABLE PARTS
INSIDE, REFER SERVICING
TO QUALIFIED PERSONNEL



FOOT
SWITCH




EXCLUSIVELY MANUFACTURED FOR		
 GIMA		GIMA SpA Via Marconi, 1 - 20060 Gessate (MI) - Italy
CUT	FREQUENCY	MANUFACTURER
200 W → 250 Ω	600 KHz	LED SpA FROSINONE - ITALY
BLEND	FREQUENCY	MODEL
120 W → 200 Ω	600 KHz	DIATERMO MB 200
FORCED COAG	FREQUENCY	CLASS
150 W → 150 Ω	600 KHz	I CF
SOFT COAG	FREQUENCY	MADE IN ITALY
90 W → 100 Ω	600 KHz	
BIPOLAR CUT	FREQUENCY	DUTY CYCLE
120 W → 250 Ω	600 KHz	int. 10 sec / 30 sec
BIPOLAR COAG	FREQUENCY	SERIAL NUMBER
100 W → 100 Ω	600 KHz	182






0051

ENSURE THAT THE CORRECT
VOLTAGE IS SELECTED BEFORE
CONNECTING TO THE MAINS
POWER SUPPLY



115 Vac - 50/60 Hz 230 Vac - 50/60 Hz
350 VA 350 VA
FUSE 2 x T 6.3 A FUSE 2 x T 3.15 A