

KIT EMERGENZA MULTIBRACCIALE - 5 bracciali
MULTI-BRACELET EMERGENCY KIT - 5 bracelets
KIT D'URGENCE MULTI-BRACELETS - 5 bracelets
KIT DE EMERGENCIA MULTIPULSERAS - 5 pulseiras
NOTFALLSET MIT MEHREREN ARMBÄNDERN - 5 Armbänder
KIT DE EMERGÊNCIA COM MÚTIPLAS PULSEIRAS - 5 pulseiras

ΚΙΤ ΕΚΤΑΚΤΗΣ ΑΝΑΓΚΗΣ ΜΕ ΠΟΛΛΑΠΛΑ ΒΡΑΧΙΟΛΙΑ - 5 βραχιόλια
Комплект за спешни случаи с няколко гривни - 5 гривни
TÖBB KARKÖTÖT SZOLGÁLÓ VÉSZHELYZETI KÉSZLET - 5 karkötő
ZESTAW AWARYJNY WIELOBRANSOLETKOWY - 5 bransoletki
KIT DE URGENȚĂ CU MULTI-BRĂȚĂRI - 5 brățări
NÖDSATS MED FLERA ARMBAND - 5 armband

CUFF SIZE - DIMENSIONE BRACCIALE

TAILLE DU BRASSARD - TAMAÑO DEL BRAZALETE - TAMANHO DO MANGUITO - MANSCHETTENGROÖßE - ΠΑΖΜΕΡ ΗΑ ΜΑΗΨΕΤΑ - ΜΕΓΕΘΟΣ ΜΑΝΣΕΤΑΣ
MANDZSETTA MÉRETE - ROZMIAR MANKIETU - DIMENSIUNEA MANȘETEI - MANSCHETTSTORLEK

XL 40.6 - 66.0 cm

COSCIA
THIGH
CUISSÉ
MUSLO
OBERSCHENKEL
COXA
ΜΗΡΟΣ
БЕДРАУН
COMB
UDO
COAPSĂ
LÅR

L 33.3 - 51.0 cm

GRANDE ADULTO
LARGE ADULT
GRAND ADULTE
ADULTO GRANDE
GROSSER ERWACHSENER
ADULTO DE GRANDE PORTE
ΜΕΓΑΛΟΣ ΕΝΗΛΙΚΟΣ
ГОЛЯМ ВЪЗРАСТЕН
NAGY FELNŐTT
DUŻY DOROSŁY
ADULT MARE
STOR VUXEN

M 25.4 - 40.6 cm

ADULTO
ADULT
ADULTE
ADULTO
ERWACHSENE
ADULTO
ΕΝΗΛΙΚΑΣ
ВЪЗРАСТЕН
FELNŐTT
DOROSŁY
ADULT
VUXEN

S 18.4 - 26.7 cm

BAMBINO
CHILD
ENFANT
NIÑO
KIND
CRIANÇA
ΠΑΙΔΙ
ДЕТЕ
GYERMEK
DZIECKO
COPIL
BARN

XS 13.0 - 19.0 cm

NEONATO
INFANT
NOURRISSON
RECIÉN NACIDO
SÄUGLING
INFANTIL
ΒΡΕΦΟΣ
КЪРМАЧЕ
CSECSEMŐ
NIEMOWLĘ
SUGARUL
SPÄDBARN

Rev.1.03.26



Gima S.p.A.
Via Marconi, 1
20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com



HONSUN (NANTONG) Co., Ltd.
N. 8, Tongxing Road
Economic&Technical Development Area
226009 Nantong City, China
Made in China

GIMA 32771

REF HS-GF-302



Shanghai International Holding Corp. GmbH (Europe)
Eiffestrasse 80
20537 Hamburg
Germany

