

# EXCLUSIVELY MANUFACTURED FOR



# GIMA

**GIMA Spa**  
Via Marconi, 1 - 20060 Gessate (MI) - Italy



**LED SpA**  
APRILIA (LT) - ITALY -

**CUT** **FREQUENCY**  
**200 W** → **250 Ω** **600 KHz**

**BLEND** **FREQUENCY**

**120 W** → **200 Ω** **600 KHz**

**ENHANCED** **FREQUENCY**

**120 W** → **250 Ω** **600 KHz**

**FORCED COAG** **FREQUENCY**

**150 W** → **150 Ω** **600 KHz**

**SOFT COAG** **FREQUENCY**

**90 W** → **100 Ω** **600 KHz**

**BIPOLAR COAG** **FREQUENCY**

**80 W** → **50 Ω** **600 KHz**

**MODEL**  
**DIATERMO MB 200 F**

**CLASS** **MADE IN**  
**I CF** **ITALY**

**DUTY CYCLE**  
**ON: 10 s / OFF: 30 s**

**SN**  
**397** \_\_\_\_\_



**CE**  
**0051**

**WARNING : NO OPERATOR SERVICEABLE PARTS INSIDE, REFER SERVICING TO QUALIFIED PERSONNEL**

ENSURE THAT THE CORRECT  
VOLTAGE IS SELECTED BEFORE  
CONNECTING TO THE MAINS  
POWER SUPPLY

115 V~ - 50/60 Hz  
350 VA - FUSE T6.3AL 250V

230 V~ - 50/60 Hz  
350 VA - FUSE T3.15AL 250V