



DIATERMO MB 122

**GIMA SpA**
Via Marconi, 1
20060 Gessate (MI)
Italy

**LED SpA**
Via Selciatella 40
04011 Aprilia (LT)
Italy

MADE IN ITALY

UDI



(01) 08033638325575
(21) 000000
(11) YYMMDD

REF GMA10100.20

SN 037 000000

 YYYY-MM

CUT
120 W => 250 Ω

BLEND
90 W => 200 Ω

FORCED COAG
80 W => 150 Ω

SOFT COAG
60 W => 100 Ω

BIPOLAR
120 W => 250 Ω

FREQUENCY
600 kHz

FREQUENCY
600 kHz

FREQUENCY
600 kHz

FREQUENCY
600 kHz

FREQUENCY
600 kHz

DUTY CYCLE
ON: 10 s / OFF: 30 s

IP32



ET037D_EN_2026-02



115V ~ - 50/60Hz
300 VA
T6.3AL 250V

