

EXCLUSIVELY MANUFACTURED FOR



GIMA

GIMA Spa
Via Marconi, 1 - 20060 Gessate (MI) - Italy



LED SpA
APRILIA (LT) - ITALY -

MODEL

DIATERMO MB 120 F

CLASS

I CF

MADE IN

ITALY

DUTY CYCLE

ON: 10 s / OFF: 30 s

SN

396 _____



CE
0051

CUT	FREQUENCY
120 W → 250 Ω	600 KHz
BLEND	FREQUENCY
90 W → 200 Ω	600 KHz
ENHANCED	FREQUENCY
90 W → 250 Ω	600 KHz
FORCED COAG	FREQUENCY
80 W → 150 Ω	600 KHz
SOFT COAG	FREQUENCY
60 W → 100 Ω	600 KHz
BIPOLAR COAG	FREQUENCY
60 W → 50 Ω	600 KHz

WARNING : NO OPERATOR SERVICEABLE PARTS INSIDE, REFER SERVICING TO QUALIFIED PERSONNEL

ENSURE THAT THE CORRECT
VOLTAGE IS SELECTED BEFORE
CONNECTING TO THE MAINS
POWER SUPPLY

115 V~ - 50/60 Hz
350 VA - FUSE T6.3AL 250V

230 V~ - 50/60 Hz
350 VA - FUSE T3.15AL 250V