

MAGNETOTERAPIA

Nome: MAG1000

Modello: NEW MAG1000

SN: 000001

Ingresso massimo: DC 4.8V, 2300mAh


Potenza d'uscita: 2x40mW

FREQUENZA portante: 27MHz

FREQUENZA modulante: 8-640Hz



IP22



I-TECH
MEDICAL DIVISION



Shenzhen Dongdixin Technology Co., Ltd.
Floor 1-2, No.3 Building, Fanshen Xusheng
Industrial Estate Xilixiaobaimang 518108
Nanshan District, Shenzhen P.R. China

MAGNETOTHERAPY

Name: MAG1000

Model: NEW MAG1000

SN: 000001

Max input: DC 4.8V, 2300mAh


Power output: 2x40mW

Carrying FREQUENCY: 27MHz

Modulating FREQUENCY: 8-640Hz



IP22



I-TECH
MEDICAL DIVISION



Shenzhen Dongdixin Technology Co., Ltd.
Floor 1-2, No.3 Building, Fanshen Xusheng
Industrial Estate Xilixiaobaimang 518108
Nanshan District, Shenzhen P.R. China

MAGNETOTERAPIA

Nombre: MAG1000

Modelo: NEW MAG1000

SN: 000001

Ingreso máximo: DC 4.8V, 2300mAh


Potencia de salida: 2x40mW

FRECUENCIA portante: 27MHz

FRECUENCIA modulante: 8-640Hz



IP22



YYYY-MM

I-TECH
MEDICAL DIVISION



Shenzhen Dongdixin Technology Co., Ltd.
Floor 1-2, No.3 Building, Fanshen Xusheng
Industrial Estate Xilixiaobaimang 518108
Nanshan District, Shenzhen P.R. China

MAGNÉTHÉRAPIE

Nom: MAG1000

Modèle: NEW MAG1000

SN: 000001

Maximale en entrée: DC 4.8V, 2300mAh


Puissance de sortie: 2x40mW

Fréquence principal: 27MHz

Modulation de fréquence: 8-640Hz



IP22



YYYY-MM

I-TECH
MEDICAL DIVISION



Shenzhen Dongdixin Technology Co., Ltd.
Floor 1-2, No.3 Building, Fanshen Xusheng
Industrial Estate Xilixiaobaimang 518108
Nanshan District, Shenzhen P.R. China

MAGNETOTHERAPIE

Name: MAG1000

Modell: NEW MAG1000

SN: 000001

Maximaler einlass: DC 4.8V, 2300mAh


Ausgangsleistung: 2x40mW

Trägerfrequenz: 27MHz

Modulationsfrequenz: 8-640Hz



IP22



I-TECH
MEDICAL DIVISION



Shenzhen Dongdixin Technology Co., Ltd.
Floor 1-2, No.3 Building, Fanshen Xusheng
Industrial Estate Xilixiaobaimang 518108
Nanshan District, Shenzhen P.R. China