

GB - Reliable digital thermometer for oral, underarm or rectal use - Fast readout - Professional accuracy - $\pm 0.1^{\circ}\text{C}$ (35.5°C-42.0°C) - $\pm 0.2^{\circ}\text{C}$ (32.0°C-35.5°C, 42.0°C-42.9°C) - Large display - Easy to use - Memory of last temperature - Storage case	DE - Zuverlässiges Digitalthermometer zum Temperaturmessen unter den Arm, in Mund oder im Rektum - Schnelles Messergebnis - Professionelle Genauigkeit - $\pm 0.1^{\circ}\text{C}$ (35.5°C-42.0°C) - $\pm 0.2^{\circ}\text{C}$ (32.0°C-35.5°C, 42.0°C-42.9°C) - Großes Display - Signalton am Ende der Messung - Speichern der letzten Temperatur - Aufbewahrungsbox	IT - Termometro digitale di precisione per uso orale, rettale ed ascellare - Misurazione rapida - Precisione elevata - $\pm 0.1^{\circ}\text{C}$ (35.5°C-42.0°C) - $\pm 0.2^{\circ}\text{C}$ (32.0°C-35.5°C, 42.0°C-42.9°C) - Ampio display - Segnale sonoro a lettura completa - Memoria ultima temperatura - Fornito con astuccio
FR - Thermomètre digital de précision à utiliser par voie orale, rectale ou axillaire. - Lecture rapide - Précision - $\pm 0.1^{\circ}\text{C}$ (35.5°C-42.0°C) - $\pm 0.2^{\circ}\text{C}$ (32.0°C-35.5°C, 42.0°C-42.9°C) - Large écran - Signal sonore après lecture - Mémoire dernière température relevée - Conditionné dans étui	ES - Termómetro digital fiable para uso oral, axilar o rectal - Lectura rápida - Precisión profesional - Exactitud elevada - $\pm 0.1^{\circ}\text{C}$ (35.5°C-42.0°C) - $\pm 0.2^{\circ}\text{C}$ (32.0°C-35.5°C, 42.0°C-42.9°C) - Display grande - Señal sonora al final de la medición - Memoriza la última temperatura - Estuche	PT - Fíavel termómetro digital para uso oral, axilar ou rectal - Leitura rápida - Grande visor - Exactidão elevada - $\pm 0.1^{\circ}\text{C}$ (35.5°C-42.0°C) - $\pm 0.2^{\circ}\text{C}$ (32.0°C-35.5°C, 42.0°C-42.9°C) - Indicador sonoro no final da medição - Memoriza a última temperatura lida - Estogo



Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com

Size 4,5x13,8x2,7 cm
S09157 - Rev. 8 - 07/21

Water-proof digital thermometer

- Termometro digitale impermeabile
- Thermomètre digital imperméable
- Digital thermometer wasserdicht
- Termómetro digital impermeable
- Termómetro digital impermeável
- Αδιάβροχο ψηφιακό θερμομέτρο
- مقياس حرارة رقمي غير نفّاذ للماء

• **Accurate for oral, rectal or underarm use**

• **Preciso per uso orale, rettale ed ascellare**

www.gimaitaly.com

REF 25562 **LOT 000013**

Gima S.p.A. - Via Marconi, 1
20060 Gessate (MI) - Italy
gima@gimaitaly.com - export@gimaitaly.com
Made in China



0476



GIMA

conf. 12 pz.
box of 12 pcs.
Boîte de 12 pcs.



0476

IP22



REF 25562

LOT 000013

TERMOMETRO DIGITALE IMPERMEABILE


W-P DIGITAL THERMOMETER °C - std. box, water-proof

W-P THERMOMÈTRE DIGITAL °C - boîte standard, étanche

TERMÓMETRO DIGITAL IMPERMEABLE °C - caja estándar, impermeab

TERMÓMETRO DIGITAL IMPERMEÁVEL °C - caixa padrão, à prova de

WASSERDICHTES DIGITALTHERMOMETER °C - Standard Packung, wass



Rev.3.07.23



Gima S.p.A. - Via Marconi, 1
20060 - Gessate (MI) - Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com
Made in China

