



ARTICOLI & APPARECCHI PER MEDICINA



Gima S.p.A - Via Marconi, 1 - 20060 Gessate (Mi) - Italy
 export@gimaitaly.com - gima@gimaitaly.com
 www.gimaitaly.com
 Made in Italy

ISTRUZIONI PER L'USO	INSTRUCTION FOR USE
• APPLICARE PRIMA DELLA SEDAZIONE • ASSICURARSI DELLA CORRETTA POSIZIONE DELL'ABBASSALINGUA NEL CAVO ORALE • APPLICARE UNA FASCIA ELASTICA PER UNA MIGLIORE STABILITÀ • NON STERILE	• INSERT BEFORE SEDATION • CHECK THAT THE TONGUE DEPRESSOR IS IN THE CORRECT POSITION • APPLY ONE HEADSTRAP FOR BETTER STABILITY • NOT STERILE
E' necessario segnalare qualsiasi incidente grave verificatosi in relazione al dispositivo medico da noi fornito al fabbricante e all'autorità competente dello Stato membro in cui si ha sede.	All serious accidents concerning the medical device supplied by us must be reported to the manufacturer and competent authority of the member state where your registered office is located.

BOCCAGLI PER ENDOSCOPIA
MOUTHPIECES FOR ENDOSCOPY

REF 71610	LOT 000009	20 PCS.	 2025-03
------------------	-------------------	---------	---



Dispositivo medico
Medical device



(01)08023279716108
 (10)000009
 (11)202503



8 023279 716108