



METRAX
GmbH Rheinwaldstr. 22
D-78628 Rottweil
Made in Germany

Rechargeable Li-Ion Battery



YYYY-MM

+EPRI738280/\$+738280000000

Type **M250/M290**



SN 738280000000

PRIMEDIC AkuPak LITE 13,2VDC / 2,5Ah / 33Wh LiFePO4

Only for use with PRIMEDIC HeartSave and
DefiMonitor XD series. Charge only with PRIMEDIC
chargers, respect charging instructions:
Max. charging voltage: 14,4VDC
Charge: 0°C ... 50°C Discharge: 0°C ... 50°C

4IXR2665-1
X=FePO4



YYYY-MM

