



1 pz.



REF 31475

LOT 0000xx

OTOSCOPIO PARKER ALOGENO - nero
PARKER HALOGEN OTOSCOPE
OTOSCOPE HALOGÈNE PARKER
OTOSCOPIO PARKER HALÓGENO

MD



GIMA S.p.A. - Via Marconi, 1 - 20060 Gessate (MI) - Italy
Made in Pakistan - www.gimaitaly.com



Rev. 1-03/2020

Size 9,8 x 2,2 x 4 cm
509102 - Rev. 5 - 03/20

Otoscopio alogeno
Parker
otoscope
Parker halogen
otoscopio halógeno Parker
otoscope halogene Parker
Parker halogen Othoskop
منظار أذن هالوجني باركر
Otosκόπιο αλογόνου Parker

Otoscopio Alogeno Parker



- Otoscope halogène Parker
- Otoscópio halógeno Parker
- Οτοσκόπιο αλογόνου Parker



- Attacco rapido a baionetta
- Lente di ingrandimento 3x mobile
- Connessione per esame pneumatico
- Illuminazione alogena 2,5V regolabile con reostato
- 3 speculum autoclavabili in dotazione Ø 2,5 - 3,5 - 4,5 mm

www.gimaitaly.com



Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com



Parker halogen otoscope



- Otoscopio halógeno Parker
- Parker halogen Othoskop
- منظار أذن هالوجني باركر



- Bayonet fitting for fast and secure attachment
- 3x movable magnifying lens
- Output for pneumatic membrane test
- 2,5V halogen light, easily adjustable
- Supplied with 3 autoclavable ear specula (Ø 2,5 - 3,5 - 4,5 mm) in a soft bag

www.gimaitaly.com



Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com

REF 31475



Gima S.p.A.
Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com
Made in Pakistan

