




ENSURE THAT THE
CORRECT VOLTAGE
IS SELECTED
BEFORE
CONNECTING TO
THE MAINS
POWER SUPPLY

MAIN POWER

230 V~ - 50/60 Hz
1000 VA

FUSE 2 X T5AL 250V

115 V~ - 50/60 Hz
1000 VA

FUSE 2 X T10AL 250V

EXCLUSIVELY MANUFACTURED FOR			
 GIMA GIMA Spa - Via Marconi, 1 20060 Gessate (MI) - Italy			
FUNCTION	MAX POWER	LOAD	FREQUENCY
PURE CUT	300 W	300 Ω	425 kHz
ENHANCED CUT	200 W	500 Ω	425 kHz
BLEND CUT	200 W	300 Ω	425 kHz
SPEEDY COAG	100 W	500 Ω	425 kHz
DEEP COAG	100 W	200 Ω	425 kHz
SPRAY COAG	70 W	2000 Ω	425 kHz
BIPOLAR CUT	95 W	150 Ω	525 kHz
BIPOLAR BLEND	90 W	150 Ω	525 kHz
BIPOLAR COAG	70 W	100 Ω	525 kHz
 LED SpA Via Salsicciello 40 APRILIA (LT) - ITALY			
MODEL			
DIATERMO MB 300 D			
DUTY CYCLE			
ON: 10 s			
OFF: 30 s			
CLASS			
I CF			
MADE IN ITALY			
SN 398 _____			

WARNING : NO OPERATOR SERVICEABLE PARTS INSIDE
REFER SERVICING TO QUALIFIED PERSONNEL

