

505.0 mm

7



STERILE FOLEY CATHETER FOR SINGLE USE

Type: 2- WAY STANDARD
Size: 12FR/CH . 5-10mL

WARNING: Do not use ointments or lubricants having a petroleum base. They will damage the catheter and cause the balloon to burst.

DIRECTIONS: After aseptic catheterization, inflate the retention balloon by inserting the syringe tip into the valve and introduce the necessary amount of sterile water.



STERILE EO



Jiangsu Suyun Medical Materials Co., Ltd.

No.1 Medicine Lane, Renmin Road, Lianyungang, Jiangsu Province, China.

Tel.:+86-518-85608151 Fax:+86-518-85466033

Made in China.

Shanghai International Holding Corp. GmbH (Europe)

Eiffestraße 80, 20537 Hamburg, GERMANY

REV.2014-10

LOT 000004

20180715

20230714