

**FOOT
SWITCH**



EXCLUSIVELY MANUFACTURED FOR



GIMA

GIMA SpA


Via Marconi, 1 - 20060 Gessate(MI)-Italy

CUT		FREQUENCY	MANUFACTURER	
160 W	→ 250 Ω	600 KHz	LED SpA	FROSINONE ITALY
BLEND		FREQUENCY	MODEL	
120 W	→ 200 Ω	600 KHz	DIATERMO MB 160	
FORCED COAG		FREQUENCY	CLASS	MADE IN ITALY
100 W	→ 150 Ω	600 KHz	I CF	
SOFT COAG		FREQUENCY	DUTY CYCLE	
80 W	→ 100 Ω	600 KHz	int. 10 sec / 30 sec	
BIPOLAR		FREQUENCY	SERIAL NUMBER	
60 W	→ 100 Ω	600 KHz	152	





ENSURE THAT THE CORRECT
VOLTAGE IS SELECTED BEFORE
CONNECTING TO THE MAINS
POWER SUPPLY



115 Vac - 50/60 Hz 230 Vac - 50/60 Hz
350 VA 350 VA
FUSE 2 x T 6.3 A FUSE 2 x T 3.15 A

WARNING : NO OPERATOR SERVICEABLE PARTS INSIDE, REFER SERVICING TO QUALIFIED PERSONNEL