




GIMA 29117

TEST CLAMIDIA - professionale
CHLAMYDIA RAPID TEST - professional
TEST DE DIAGNOSTIC DE CHLAMYDIA
TEST CLAMIDIA - profesional



6 936983 186971

Importato da/Imported by/Importé par/Importado por:
GIMA S.p.A. - Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com



29117
000027

SCATOLA
BOX



RACCOLTA CARTA
PAPER COLLECTION

BUSTINA
POUCH



RACCOLTA
INDIFFERENZIATA
MIXED WASTE COLLECTION

ISTRUZIONI
PER L'USO
USER INSTRUCTIONS



RACCOLTA CARTA
PAPER COLLECTION

Verifica le disposizioni del tuo comune per la gestione dei rifiuti
Please refer to the provisions of your municipality
concerning waste disposal

Chlamydia Rapid Test Cassette (Swab/Urine)/ Cassetta per test rapido per Clamidia (Tampone/Urina)

Contents/Contenuto:

Test Cassettes/Cassetta del Test - 20

Extraction Tubes/Provetta di estrazione - 20

Sterile Swabs/Tamponi sterili - 20

Extraction Reagent 1/Reagente di estrazione 1 - 1

Extraction Reagent 2/Reagente di estrazione 2 - 1

Package Insert/Foglio illustrativo - 1

Dropper Tips/Puntas contagocce - 20

Workstation/Postazione di lavoro - 1

REF ICH-502


LOT



CE 0123



IVD



EC REP

MedNet EC-REP GmbH

Borkstrasse 10

48163 Muenster Germany



Hangzhou AllTest Biotech Co., Ltd.

#550, Yinhai Street,

Hangzhou Economic & Technological Development Area,

Hangzhou, 310018 P.R. China

Web: www.alltests.com.cn

Email: info@alltests.com.cn

MADE IN CHINA