

Blood Lancets

- Blodlancetter
- Lancete za kri
- Lancety na odber krvi
- Ланцети за кръв
- Kraujo lancetai
- Asins lancetes
- Verevõtulantsetid
- Ланцеты для кровопускания
- واخزات الدم




28G

32G

30G

Instructions for use - Istruzioni per l'uso :



Twist cap
Ruotare
il tappo

Pull cap
Tirare
il tappo

100 pezzi
pcs.

Information about the matching lancing device
Informazioni sui dispositivi pungidito compatibili

- LINKFAR Lancing Devices
- Lianfa Lancing Devices
- OneTouch® UltraSoft®
- OneTouch® UltraVue
- FreeStyle Lite®
- Bayer Microlet®
- Bayer Microlet® Vaculance
- Glucolet Lancing Device
- BD Lancet Device
- Easy Touch® lancing device
- Mylife Unio (Ypsomed)
- Mylife Pura (Ypsomed)
- Beurer Lancing Device
- Autolet Lancing Device
- Glucolet Lancing Device
- SD Biosensor Lancing Device
- TaiDoc Lancing Device
- CareSoft Lancing Device
- Tempo Lancing Device
- Autolet
- MICROLET® NEXT
- Medifun
- SteriLance
- AIMSCO Lancing device

Lancette pungidito

- Lancettes
- Lancetas de sangre
- Lancetas para coleta de sangue
- Lancetas para coleta de sangue
- Lancety do pobierania próbek krwi
- Σκαρφιστήρες
- Lancety pro odběr krve
- Blodlansetter
- Verilansetit
- Lanžete recoltare sânge
- Lancete za krv
- Ujjbegyszűrő lándzsa
- Blodlansetter
- Bloedprikpennen




28G

32G

30G

100 pezzi
pcs.

- ☐ 28G **GIMA 23916**
- ☐ 30G **GIMA 23918**
- ☐ 32G **GIMA 23921**



SIZE: 75*50*92mm

Shandong Lianfa Medical Plastic Products Co., Ltd.
No.1 Shuangshan Sanjian Road, Zhangqiu,
Jinan City, 250200, Shandong
P.R. China - Made in China

REF 01-1328-II, 01-1330-II, 01-1332-II

EU REP Linkfar Healthcare GmbH
Niederrheinstraße 71, 40474 Düsseldorf, Germany
UK REP MedPath Limited, 27 Old Gloucester Street,
London, WC1N 3AX United Kingdom

Gima S.p.A.
Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com

CE 0197



Open from the other side

S09691 Rev. 2 - 05/26
Size 7,5x5x9,2 cm