



Jiangsu Suyun Medical Materials Co., Ltd  
江苏苏云医疗器材有限公司

Add: No.1, Medicine Lane, Renmin Road, Lianyungang, Jiangsu, China  
地址: 中国江苏省连云港市人民路制药巷1号 (222000)

Tel: +86-518-85608151/85476682

Fax: +86-518-85451119/85450398

www.suyunmedical.com

## Declaration of DEHP FREE

To whom it may concern,

DEHP FREE is a concept, a limit, not 0% in the product. Europe, US, China, as well as other countries like Brazil, etc. have the related standard, which stipulate the DEHP content should not more than 1000ppm in the product, means its proportion should not more than 0.1%.

We have authorized the Jiangsu Province Institute for Medical Equipment testing to calculate the DEHP content in our raw material. The test of this institute is acceptable by Chinese SFDA. Test method is according to EN14372 item 6.3.2. They put 2g raw material into 50ml reagent, final result is 0.21mg/100mL, means 0.105mg/50mL. In the 2g plastic granula, 0.105mg takes only 52.5ppm, far less than the 1000ppm of the standard.

So we can say that the DEHP FREE content of our following PVC products is acceptable according to the requirement of international standards.

Sterile Nelaton Catheters code. F5110 F5100  
Sterile Feeding Tubes code. F5400  
Sterile Stomach tubes code. F5452 F5411  
Sterile Suction Catheters code. F5200 F5220 F5202

Jiangsu Suyun Medical Materials Co., Ltd  
May. 13rd 2016

