

SIMPLE WHITE SMALL - EMERGENCY TROLLEY IN COATED STEEL WITH 5 DRAWERS AND ACCESSORIES

Cod. KSM-EM1S

TECHNICAL DATA SHEET:

CONSTRUCTION: Solid, sturdy, made of powder epoxy coated steel and techno-polymer (plastic material) .
Maximum lightness and maneuverability.

DESIGN: Original design and great attention to every little detail. Suitable for any hospital department. The special design makes trolley cleaning easy.

STRUCTURE: Self-supporting and monocoque.

TOP SURFACE: Comfortable large working surface with small objects compartment, made of techno-polymer with no sharp edges, in order to provide maximum safety and cleanability. Height from floor: mm 1007

PUSH HANDLE: ergonomic.

FRONT SIDE: Pre-settled for central locking system and modular drawers, without sharp edges, entirely made of techno-polymer. Silent running sliding drawers with bearings and security lock for central locking.

DEFIBRILLATOR HOLDER: Plexiglass tray (dim. mm 390 x 320 x 20 H) and stainless steel tubular support. Able to rotate 360° with possibility to lock every 10° and equipped with 2 safety straps. Straps locking device suitable for any unit size.

I.V. POLE: Made of stainless steel, adjustable height.

CPR MASSAGE BOARD : Made of transparent plexiglass.

DRAWERS: Made of Powder epoxy powder steel, n°2 dimension drawer mm. 450x475x75H, n°2 dimension drawer mm. 450x475x151H and n°1 dimension drawer mm. 450x475x228H. Colored dot stickers are included for an easy identification.

LOCKING SYSTEM: Safe locked drawers by the CPR board and plastic seals (code K617363), single use. N. 20 spare plastic seals are included and provided with the trolley.

LEFT SIDE : additional plastic working table; built-in compartment with 4 removable shelves, made of plastic techno polymer.

RIGHT SIDE: Built-in compartment made of techno polymer with n°1 large waste container 15lt and n°1 catheter holder.

BACKSIDE : oxygen tank holder compartment with fixing straps, n° 1 holding rail with hooks.

CASTOR: Set of 4 swivel castors (diam.125 mm), integrated thread guard, which n° 2 conductive wheels with lever locking system, double ball bearing.

BUMPER: integrated base bumper and 4-corner bumper protection on each castor.

SAFETY: in accordance with the current international regulations.

USE: Ideal for any emergency and first aid.

DIMENSIONS: mm. 850 x 714 x 1182H

