



SavePads

Multifunction electrodes for adults



KEY FEATURES

- Disposable, self-adhesive multifunction electrodes for adults for defibrillation, pacing, monitoring and cardioversion
- 100 % quality control
- For up to 50 shocks (360 J biphasic)
- For up to 24 hours of monitoring
- For up to 1 hour of pacing using 140mA / 120bpm (pulse duration 20ms)
- For up to 8 hours of pacing using 70mA / 60bpm (pulse duration 20ms)
- Inspection of pads every 30 minutes

PRODUCT

Characteristics: disposable, pre-gelled, self-adhesive, single use, non-sterile
Size: adult
Shelf life: 36 months unopened pouch
Temperature conditions
Transport, max. 10 days: -40°C - 75°C
Storage: 0°C - 50°C
Operating: 0°C - 50°C
Humidity (no condensation): 5% - 95%
Atmospheric Pressure: 500hPa - 1.060hPa

CLASSIFICATION

Classification (MDD Annex IX): Class I

PAD DIMENSIONS

Electrode Shape: rectangular
Total Area: appr. 148 cm²
Adhesive Area: appr. 145 cm²
Gel Area / Active Area: appr. 87 cm²
Gel Thickness: 0.60 ± 0.10 mm

MATERIALS PAD

Backing Material: Adhesive PE Foam
Conductive Material: Tin
Gel: Adhesive hydrogel
Release Liner: Siliconized PET foil

ELECTRODE PACKAGING

Pieces / Pouch: 1 Set (2 pads / pouch)

MATERIALS PACKAGING

Pouch: PET / Al / PE (PET Polyester, Al Aluminium, PE Polyethylene)

REACH CANDIDATE LIST SUBSTANCES

Information according to Art. 33 for substances contained above 0.1 wt%

Substance: Boric acid*

*Based on the information currently available, this substance does not pose any risk if the article is used as intended. Additional information is not necessary to allow safe use of the article.

CONNECTION

Connector type:



Cable length: 1,2 m

ACCESSORIES

SavePads AED: Resuscitation kit: Artificial respiration cloth, razor, disposable gloves

Metrax GmbH
Rheinwaldstr. 22
D-78628 Rottweil
Germany

info@primedic.com
www.primedic.com

Release date: 02/2019
24034 | EN | A

www.primedic.com



PRIMEDIC
Saves Life. Everywhere.