



SavePads PreConnect

Multifunction electrodes for adults / pediatrics



KEY FEATURES

- Disposable, self-adhesive multifunction electrodes for adults and pediatrics for defibrillation, pacing, monitoring and cardioversion
- 100 % quality control
- For up to 50 shocks (360 J biphasic)
- For up to 24 hours of monitoring
- For up to 1 hour of pacing using 140mA / 120bpm (pulse duration 20ms)
- For up to 8 hours of pacing using 70mA / 60bpm (pulse duration 20ms)
- Inspection of pads every 30 minutes

PRODUCT

Characteristics: disposable, pre-gelled, self-adhesive, single use, non-sterile
Size: adult / child
Shelf life: 36 months unopened pouch
Temperature conditions
Transport, max. 10 days: -40°C - 75°C
Storage: 0°C - 50°C
Operating: 0°C - 50°C
Humidity (no condensation): 5% - 95%
Atmospheric Pressure: 500hPa - 1.060hPa

CLASSIFICATION

Classification (MDD Annex IX):
 Class I

PAD DIMENSIONS

Electrode Shape: rectangular
Total Area: appr. 148 cm²
Adhesive Area: appr. 145 cm²
Gel Area / Active Area: appr. 87 cm²
Gel Thickness: 0.60 ± 0.10 mm

MATERIALS PAD

Backing Material: Adhesive PE Foam
Conductive Material: Tin
Gel: Adhesive hydrogel
Release Liner: Siliconized PET foil

ELECTRODE PACKAGING

Pieces / Pouch: 1 Set (2 pads / pouch)

MATERIALS PACKAGING

Pouch: PET / Al / PE (PET Polyester, Al Aluminium, PE Polyethylene)

CONNECTION

Connector type: pre-connectable



Cable length: 2 m

ACCESSORIES

SavePads PreConnect AED / AS:
 Resuscitation kit: Artificial respiration cloth, razor, disposable gloves, scissors

Metrax GmbH
 Rheinwaldstr. 22
 D-78628 Rottweil
 Germany

info@primedic.com
 www.primedic.com

Release date: 02/2019
 24033 | EN | A

www.primedic.com



PRIMEDIC
 Saves Life. Everywhere.