

AESCULAP TC SURGICAL INSTRUMENTS - DUROTIP® AND DUROGRIP®

SCISSORS

39351

- **39350** Mayo TC durotip thin scissors - straight - blunt/blunt - 17 cm - BC242R
- **39351** Mayo TC durotip thin scissors - curved - blunt/blunt - 17 cm - BC243R
- **39355** Mayo TC durotip scissors - straight - blunt/blunt - 17 cm - BC252R
- **39356** Mayo TC durotip scissors - curved - blunt/blunt - 17 cm - BC253R

39366

- **39360** Baby metzenbaum TC durotip scissors - curved - blunt/blunt - 14.5 cm - BC259R
- **39365** Metzenbaum TC durotip scissors - curved - blunt/blunt - 14.5 cm - BC261R
- **39366** Metzenbaum TC durotip scissors - curved - blunt/blunt - 18 cm - BC263R
- **39367** Metzenbaum TC durotip thin scissors - curved - blunt/blunt - 18 cm - BC271R

NEEDLE HOLDER

39391

- **39390** Castroviejo TC durogrip needle holder - straight - 0.2/14.5 cm - BM002R
- **39391** Castroviejo TC durogrip needle holder - straight - 0.4/18 cm - BM003R
- **39392** Castroviejo TC durogrip needle holder - straight - 0.2/18 cm - BM004R

39396

- **39395** Baby-crile-wood TC durogrip needle holder - straight - 15 cm - BM013R
- **39396** Baby-crile-wood TC durogrip needle holder - straight - 14.5 cm - BM016R
- **39400** Crile-wood TC durogrip needle holder - straight - 18.5 cm - BM017R

AESCULAP® - a B. Braun brand

B | BRAUN
SHARING EXPERTISE

39406

- **39405** Hegar-mayo TC durogrip needle holder - straight - 15 cm - BM065R
- **39406** Hegar-mayo TC durogrip needle holder - straight - 18.5 cm - BM066R

39411

- **39410** Hegar-olsen TC durogrip needle holder - straight - 16.5 cm - BM124R
- **39411** Hegar-olsen TC durogrip needle holder - straight - 14.5 cm - BM128R

39415

- **39415** Mathieu tc durogrip needle holder - straight - 17.3 cm - BM151R

QUICK TOURNIQUETS

25748

Multilingual box: GB, FR, IT, ES, PT, DE, NL, DK, PL, SE, CZ, HU, GR, RO, HR



ABS buckle



One hand use and quick release tourniquet

GIMA code	TOURNIQUETS	Material	Colour	Length cm	Width cm	Fastening systems	Origin
25748	Quick tourniquet	Spandex/Elastan	Blue	45	2.5	One hand use - quick release	Asia
25749	Quick tourniquet - latex free	Spandex/Elastan	Green	45	2.5	One hand use - quick release	Asia

CUFFED AND UNCUFFED TRACHEOSTOMY TUBES

CUFFED AND UNCUFFED TRACHEOSTOMY TUBES

Curve tubes designed to rapidly provide an airway within seconds during an emergency by inserting them into a tracheostomy stoma, a hole created in the neck and windpipe.

The hollow tube whether with or without cuff, is directly placed into the trachea through a surgical incision, using a wire-guided progressive dilation technique.

- tube is made of medical grade, non-toxic, clear and kink-resistant PVC to protect delicate mucosal tissue
- smooth rounded obturator's tip reduces trauma during insertion

- neck plate with printed size for easy identification
- standard connector

Tube includes: cuff (cuffed only), inflation line, valve, pilot balloon (cuffed only), obturator, neck plate, connector and tube ties.

Sterile, single patient.

Class IIa. DEHP-free.

Made with no natural rubber or latex.

Cuffed tracheostomy tubes additional features

- high-volume low-pressure cuff provides an effective low-pressure seal
- soft, thin-walled cuff ensures effective sealing and atraumatic intubation and extubation
- efficient one-way valve for cuff inflation and deflation



57671-3

Tracheostomy tube without cuff

57679-57685

Tracheostomy tube with cuff

All tracheostomy tubes are individually packed in blister.

GIMA code	UNCUFFED TRACHEOSTOMY TUBES	Minimum order
57671	Uncuffed tracheostomy tubes Ø 4.0 mm	box of 10
57672	Uncuffed tracheostomy tubes Ø 4.5 mm	box of 10
57673	Uncuffed tracheostomy tubes Ø 5.0 mm	box of 10
CUFFED TRACHEOSTOMY TUBES		
57679	Cuffed tracheostomy tubes Ø 6.0 mm	box of 10
57680	Cuffed tracheostomy tubes Ø 6.5 mm	box of 10
57681	Cuffed tracheostomy tubes Ø 7.0 mm	box of 10
57682	Cuffed tracheostomy tubes Ø 7.5 mm	box of 10
57683	Cuffed tracheostomy tubes Ø 8.0 mm	box of 10
57684	Cuffed tracheostomy tubes Ø 8.5 mm	box of 10
57685	Cuffed tracheostomy tubes Ø 9.0 mm	box of 10