

S/S AISI 202 "GOLDEN" KIDNEY DISHES

Golden kidney dishes easy to clean and sterilize. Nickel and chromium free, makes it hypoallergenic. Coating prevents corrosion.

GIMA code	GOLDEN S/S KIDNEY DISHES	Capacity
37750	Tray 6" - 150 mm	200 ml
37751	Tray 8" - 200 mm	500 ml
37752	Tray 10" - 250 mm	800 ml
37753	Tray 12" - 300 mm	1500 ml



S/S AISI 304 BOWL-TAPERS S/S AISI 304 BEAKERS

Autoclavable at 121° C



GIMA code	MIXING BOWL-TAPER	Size mm	Capacity
37800	Mixing bowl-taper	Ø80xh 33	100 ml
37802	Mixing bowl-taper	Ø105xh 53	300 ml
37804	Mixing bowl-taper	Ø285xh 120	4000 ml

Autoclavable at 121° C
Thickness 0.6 mm



GIMA code	BEAKERS	Size mm	Capacity
37810	Beaker	Ø70xh 82	250 ml
37812	Beaker with handle	Ø105xh 130	1000 ml
37813	Beaker with handle	Ø130xh 168	2000 ml

DEEP KIDNEY DISHES



GIMA code	AISI 202 DEEP KIDNEY DISHES	Capacity
37865	Kidney dish - 200 mm	500 ml
37866	Kidney dish - 250 mm	800 ml
37867	Kidney dish - 300 mm	1500 ml

SHALLOW KIDNEY DISH



GIMA code	AISI 202 SHALLOW KIDNEY DISHES	Capacity
37860	Kidney dish - 150 mm	200 ml
37861	Kidney dish - 200 mm	500 ml
37862	Kidney dish - 250 mm	800 ml
37863	Kidney dish - 300 mm	1500 ml

GALLIPOTS



GIMA code	AISI 202 GALLIPOTS WITH SPOUT	Capacity
37851	Gallipot with spout - 60x30 mm	60 ml
37852	Gallipot with spout - 80x40 mm	100 ml
37853	Gallipot with spout - 100x50 mm	250 ml
37855	Gallipot with spout - 150x65 mm	700 ml

INSTRUMENT TRAYS



GIMA code	AISI 202 INSTRUMENT TRAYS
37870	Instrument tray - 200x150xh 40 mm
37871	Instrument tray - 220x150xh 30 mm
37872	Instrument tray - 220x150xh 50 mm
37874	Instrument tray - 250x200xh 50 mm
37875	Instrument tray - 300x210xh 40 mm
37877	Instrument tray - 350x260xh 40 mm

S/S STUDENT DISSECTION KIT

• 26955 S/S AISI 410 STUDENT DISSECTION KIT - 11 pieces
Standard dissecting kit suitable for medical students.
Kit includes 11 instruments and case.

KIT CONTENT	
Scalpel knife handle no. 3	Teasing needle straight plastic handle
Iris scissors 4.5"	Teasing needle curved plastic handle
Mosquito forceps 5"	Dissecting forceps curved
Tissue forceps tooth 1x2	Operating scissors Sh/Bl straight 5"
Mall probe	Dressing forceps
Groove director 5.5"	

