

BD INSYTE™ AUTOGUARD™ SHIELDED IV CATHETERS

BD Insyte™ Autoguard™ shielded IV catheters are over-the-needle, intravascular (IV) catheters that incorporate advanced safety features including a push button that instantly retracts the needle.

Enhanced Safety

BD push-button needle shielding technology is demonstrated to reduce risk of accidental needlestick injuries.

Insertion success in sight

BD Instaflash™ needle technology incorporates a notched needle, which clinically demonstrated improved first-stick insertion success, reducing painful hit-and-miss insertions.

Reduced complications

BD Vialon™ biomaterial softens up to 70% in the vein, enabling longer dwell times and reducing the chance of mechanical phlebitis.

Intended to be inserted into a patient's peripheral vascular system for short term use to sample blood, monitor blood pressure, or administer fluids. Available with or without wings. Multilingual manual and box in GB, FR, IT, ES, DE, PT, NL, DK, SE, NO, FI, LT, EE, PL, CZ, SK, HU, SI, HR, RO, BG, GR, RU, TR, Arabic.



23902



23901

GIMA code	BD code	BD INSYTE™ AUTOGUARD™			Minimum order
		Gauge	Gauge mm	Colour	
23900	381844	18G	1.3x30 mm	green	box of 50
23901	381834	20G	1.1x30 mm	pink	box of 50
23902	381823	22G	0.9x25 mm	blue	box of 50
23903	381812	24G	0.7x19 mm	yellow	box of 50
BD INSYTE™ AUTOGUARD™ - WINGED					
23905	381944	18G	1.3x30 mm	green	box of 50
23906	381934	20G	1.1x30 mm	pink	box of 50
23907	381923	22G	0.9x25 mm	blue	box of 50
23908	381912	24G	0.7x19 mm	yellow	box of 50

BD INSYTE™ IV CATHETERS

BD Insyte™ IV catheters are intravascular radiopaque catheters with Luer-Lock cone.

Intended to be inserted into patient's peripheral vascular system for short term use to sample blood, monitor blood pressure or administer fluids. Needle is internally and externally lubricated with triple sharpening (backcut grinding). BD Vialon™ material softens up to 70% in the vein, enabling longer dwell times and reducing the chance of mechanical phlebitis and infiltration. Class IIa.

Multilingual manual and box: GB, FR, IT, ES, PT, DE, NL, DK, SE, NO, FI, IS, LT, LV, EE, PL, CZ, SK, HU, SI, HR, RO, BG, GR, UA, TR, RU, Arabic, KZ, PK, CN, KR.



23888

GIMA code	BD code	BD INSYTE™			Minimum order
		Gauge	Gauge mm	Colour	
23885	381244	18G	1.31x30 mm	green	box of 50
23887*	381234	20G	1.09x30 mm	pink	box of 50
23888*	381223	22G	0.88x25 mm	blue	box of 50
23889*	381212	24G	0.7x19 mm	yellow	box of 50

*BD Instaflash™ needle technology incorporates a notched needle, which clinically demonstrated improved first-stick insertion success, reducing painful hit-and-miss insertions.

TERUMO® SURFLO™ INFUSION SETS

Surflo™ Winged Infusion Sets (Scalp Vein or SV set) are intended to introduce fluids into the vascular system after venipuncture. To ensure stability and grip during venipuncture, they are equipped with color-coded, interlocking wings. Thanks to the slightly angled wings (<3° in relation to the cannula), less pressure is put on the blood vessel when puncturing. They allow smooth skin penetration. Not made with DEHP and have no natural rubber latex components. Tube length: 30 cm. Class IIa, Medical Device.

Multilingual box: GB, FR, IT, ES, PT, DE, PL, NL, HU, GR, BG, RO, FI, NO, SE, DK, CZ, SK, SI, RS, EE, LV, LT, TR, RU.



37594



SURFLO™ WINGED INFUSION SETS - STERILE				
GIMA code	TERUMO® code	Needle	Colour	Minimum order
37591	SV*19NL30	19G - 1.1x19 mm	ivory	box of 50
37592	SV*21NL30	21G - 0.8x19 mm	green	box of 50
37593	SV*22NL30	22G - 0.7x19 mm	black	box of 50
37594	SV*23NL30	23G - 0.6x19 mm	blue	box of 50
37595	SV*25NL30	25G - 0.5x19 mm	orange	box of 50

TERUMO® SURFLO™ I.V. CATHETERS

Surflo™ I.V. Catheters are designed to ensure a stable form (lower risk of split tips and peel back) during punctures under difficult conditions. They are made of medical grade EFTE. Class IIa, Medical Device.

Multilingual box: GB, FR, IT, ES, PT, DE, PL, NL, HU, GR, BG, RO, FI, NO, SE, DK, CZ, SK, SI, RS, EE, LV, LT, TR, RU.



37587



SURFLO™ I.V. CATHETERS - STRAIGHT - STERILE				
GIMA code	TERUMO® code	Needle mm	Colour	Minimum order
37580	SR+OX1464C1	14Gx2 1/2" - 2.15x64	orange	box of 50
37584	SR+OX1864C1	18Gx2 1/2" - 1.3x64	green	box of 50
37585	SR+OX2032C1	20Gx1 1/4" - 1.1x32	pink	box of 50
37587	SR+OX2225C1	22Gx1" - 0.85x25	blue	box of 50
37588	SR+OX2419C1	24Gx3/4" - 0.67x19	yellow	box of 50

TERUMO® SURSHIELD™ VERSATUS I.V. CATHETERS WITH PROTECTOR

Surshield™ Versatus I.V. Catheter with protector allows to completely close the needle, in order to prevent and blood contamination and accidental needle stick injuries, when disposing of the intravenous cannula. Sharps protection

SURSHIELD™ VERSATUS I.V. CATHETERS - STERILE				
GIMA code	TERUMO® code	Needle	Colour	Minimum order
37571	SR+DS1645PD	16G - 1.75x45 mm	grey	box of 50
37573	SR+DS1832PD	18G - 1.3x32 mm	green	box of 50
37574	SR+DS1845PD	18G - 1.3x45 mm	green	box of 50
37575	SR+DS2032PD	20G - 1x32 mm	pink	box of 50
37576	SR+DS2225PD	22G - 0.9x25 mm	blue	box of 50
37577	SR+DS2419PD	24G - 0.74x19 mm	yellow	box of 50
37578	SR+DS2619PD	26G - 0.64x19 mm	violet	box of 50

for healthcare workers provides an important tool in the healthcare setting. Special tip design for smooth and less painful vein-puncture. Class IIa, medical device.

Multilingual box: GB, FR, IT, ES, PT, DE, PL, NL, HU, GR, BG, RO, FI, NO, SE, DK, CZ, SK, SI, RS, EE, LV, LT, TR, RU.



37575

