

BD VENFLON™ PRO SAFETY I.V. CATHETERS - STERILE

BD Venflon™ Pro Safety Needle Protected IV Cannula is a ported, passive safety peripheral IV catheter (PIVC) designed to help minimise the risk of unanticipated blood spatter and needlestick injuries (NSIs).

BD Vialon™ catheter material (polyurethane, radiopaque and biocompatible) softens in the vein up to 70%, enabling longer dwell times, minimizing restarts and reducing the risk of mechanical phlebitis. Ergonomic design to accommodate a variety of techniques.


Back-cut needle design to offer optimal skin and vein penetration.

Flexible wings help provide comfortable securement for the patient. **"Snap" cap** helps prevent port being accidentally opened.

Removable white vent cap provides flexibility to address a variety of needs.

GIMA code	BD code	BD VENFLON PRO SAFETY NEEDLES			Minimum order
		Gauge	Gauge mm	Colour	
23712	393226	18G	1,3x32 mm	green	box of 50
23713	393224	20G	1,1x32 mm	pink	box of 50
23714	393222	22G	0,9x25 mm	blue	box of 50

Multilingual instructions in:
GB, FR, IT, ES, PT, DE, NL, PL,
DK, FI, SE, NO, EE, LV, LT, CZ,
SI, SK, GR, BG, HR, HU, RO,
UA, TR, RU, KZ, Arabic.



BD CONNECTA™ PLUS 3

3-way stopcocks for intravenous therapy and pressure monitoring, with rotating Luer lock. Single use, sterile. **CE 2797**


- high flow capacity - full 360° rotation
- low priming volume (dead space): 0.23 ml
- lipid resistant: material allows the use of lipids and chemotherapeutic agents



Multilingual manual: GB, FR, IT, ES, PT, DE, NL, DK, SE, NO, FI, EE, PL, CZ, SK, HU, SI, HR, RO, BG, GR, UA, TR, RU, KR

GIMA code	BD code	BD CONNECTA™ PLUS 3 THREE WAY STOPCOCKS	Minimum order
23854	394600	3-way stopcocks	box of 100
23855	394995	3-way stopcocks with 10 cm extens. cable	box of 50
23856	394951	3-way stopcocks with 50 cm extens. cable	box of 20

STERILE SURGERY SETS



- **23639 UNIVERSAL SURGICAL SET - sterile - box of 12**
- 1 table cover for instruments with absorbent band 160x200 cm
- 1 Mayo table cover with absorbent band 80x145 cm
- 1 standard surgical drape with adhesive tape 150x240 cm
- 2 standard surgical drapes with adhesive tape 75x90 cm
- 1 standard surgical drape with adhesive tape 150x200 cm
- 4 towels 30x40 cm
- 2 adhesive strips 8x50 cm

BI-LAYER NON-WOVEN STERILE SURGICAL DRAPES

Single-use sterile surgical drapes made of non-woven, absorbent and waterproof material. Made with no natural rubber or latex. Class IIa. **MADE IN ITALY**



Multilingual instructions and label: GB, FR, IT, ES, PT, DE, NL, DK, SE, NO, FI, LT, LV, EE, PL, CZ, SK, HU, SI, HR, RO, BG, GR, Arabic

GIMA code	DRAPES	Minimum order
23580	Drape 50x50	box of 350
23581	Drape 45x75	box of 320
23582	Drape 75x90	box of 170
23583	Drape 100x100	box of 120
23584	Drape 100x150	box of 80
23585	Drape 150x150	box of 70
23586	Drape 150x200	box of 50
23587	Drape 150x240	box of 40
23588	Drape 200x200	box of 40
23589	Drape 200x250	box of 30


LEGGINGS

- **23658 PAIR OF LEGGINGS 79x129 cm - sterile - box of 40**

For urology and gynaecology.

MADE IN ITALY

Multilingual instructions and label:
GB, FR, IT, ES, PT, DE, NL, DK, SE,
NO, FI, LT, LV, EE, PL, CZ, SK, HU,
SI, HR, RO, BG, GR, Arabic



STOCKINETTE

- **23659 STOCKINETTE LEGGING 25x75 cm - medium - sterile - box of 60**
- **23660 STOCKINETTE LEGGING 35x120 cm - large - sterile - box of 40**

For orthopedics and gynaecology.

Multilingual instructions and label: GB, FR, IT, ES, PT, DE, NL, DK, SE, NO, FI, LT, LV, EE, PL, CZ, SK, HU, SI, HR, RO, BG, GR, Arabic

MADE IN ITALY



SURGERY TREATMENT SET - PROFESSIONAL

- **26952 SURGERY TREATMENT SET - sterile - box of 30**

Single-use sterile set intended for use by qualified healthcare professionals for suturing, tissue clamping, skin incision, and other procedures. For use in hospitals and medical practices. The instruments provide the same quality as reusable ones. Medical device, Class IIa. Multilingual manual GB, FR, IT, ES, DE.

KIT CONTENT

- 1 scalpel n. 15 Swann-Morton
- 1 Micro Halsey needle holder 13 cm
- 1 Adson forceps ultra thin 1x2 teeth
- 1 curved blunt iridectomy forceps
- 1 absorbent laminated patient drape with holes 50x60 cm, Ø 8 cm
- 1 absorbent table drape 38x45 cm
- 5 non-woven swabs 7.5x7.5 cm, 30 g
- 1 three-piece syringe 2 ml
- 1 hypodermic needle 25G x 16 mm
- 1 hypodermic needle 21G 1/2 x 40 mm

