

STERI-STRIP™ - REINFORCED ADHESIVE SKIN CLOSURES

Reinforced adhesive skin closures are porous non-woven adhesive strips indicated for use as a skin closure device in the treatment of lacerations and surgical incisions. Designed to let the skin breathe and reinforced with filaments for strength.

- comfortable to wear compared to sutures

- non-invasive and sterile design

helps to reduce risk of infections

Pressure sensitive hypoallergenic adhesive.

Multilingual instructions: GB, FR, IT, ES, PT, DE, NL, DK, SE, NO, FI, GR.

Box: FR, IT, DE.

GIMA code	3M™ code	STERI-STRIP™ - STERILE	Minimum order
25270	1540NP	3x75 mm - white	box of 2
25271	1541NP	6x75 mm - white	box of 2
25272	1546NP	6x100 mm - white	box of 1
25273	1547NP	12x100 mm - white	box of 1
25275	1551NP	6x75 mm - skin	box of 2



25275 Skin Steri-Strip



25270 White Steri-Strip

MEDIPORE™ AND TEGADERM™ - WOUND CARE

Non-adherent absorbent pad bonded to a larger soft cloth backing coated with a border of hypoallergenic, water-resistant adhesive.

Mechanical barrier for absorption of exudates for acute wounds

such as cuts, burns,

abrasions, IV catheter sites and surgical

incisions. Sterile, highly breathable and conformable. It does not contain latex.

Multilingual instructions: GB, FR, IT, ES, PT,

DE, NL, DK, SE, NO, FI, LT, LV, EE, PL, CZ,

SK, HU, RS, SI, HR, RO, BG, GR, UA, RU, TR.



25277

Transparent sterile dressing with absorbent, non-stick pad. Breathable and waterproof polyurethane coating ensures total protection from external contamination

and reduces dressing

changes, allowing patients

to bathe and shower. High absorption

capacity and optimal skin transpiration,

accelerating wound healing.

Multilingual instructions: GB, FR, IT, ES, PT,

DE, NL, DK, SE, NO, FI, GR.



25282

Transparent film dressing that provides an effective protective barrier against external contaminants.

Ideal for treating

superficial wounds, first

and second degree burns

and for covering and fixing

devices to the skin. It

protects from dirt, bacteria

and viruses. Practical side

wings for easy application.

Waterproof and latex free.

Multilingual instructions: GB, FR, IT, ES,

PT, DE, NL, DK, SE, NO, FI, LT, LV, EE, PL,

CZ, SK, HU, RS, SI, HR, RO, BG, GR, UA,

TR, RU.



25285

GIMA code	3M™ code	MEDIPORE™ + PAD - STERILE	Minimum order
25277	3562NP	5x7.2 cm	box of 5
25278	3566NP	10x10 cm	box of 5
25279	3569NP	10x15 cm	box of 5
25280	3570NP	10x20 cm	box of 5

GIMA code	3M™ code	TEGADERM™ + PAD - STERILE	Minimum order
25281	3582NP	5x7 cm	box of 5
25282	3586NP	9x10 cm	box of 5
25283	3584NP	6x10 cm	box of 5
25284	3589NP	9x15 cm	box of 5

GIMA code	3M™ code	TEGADERM™ FILM - STERILE	Minimum order
25285	1624NP	6x7 cm	box of 5
25286	1626NP	10x12 cm	box of 5

MICROPORE™ AND TRANSPORE™ - TAPES FOR MEDICAL FIXATION

Tapes with silicone adhesive, particularly gentle on the skin that guarantees a painless removal.

For sensitive skin.

It stays on the skin for up to 48 hours.

Bi-directional tearing without the use of scissors.

Suitable for infants and elders.

Multilingual box

and instructions:

IT, FR, DE.



25289

Highly breathable and gentle on the skin tapes that offer a multifunctional fixation, ideal for fixing dressings. On the back of the support there's a substance called LAB (Low Adherence

Backing) for

low adhesion

to the back.

Hypoallergenic

adhesive. Box,

instructions:

IT, FR, DE.



25291

Transparent, porous and breathable surgical tape with optimal adhesion to skin and medical devices, ideal for fixing dressings, tubes and other devices.

- easy to tear

- easy to use with gloves

- waterproof

- hypoallergenic and latex-free

Multilingual box and instructions:

IT, FR, DE.



25293

GIMA code	3M™ code	MICROPORE™ S - STERILE	Minimum order
25289	2775NP-1	2.5 cm x 5 m	box of 1

GIMA code	3M™ code	MICROPORE™ - STERILE	Minimum order
25287	1533NP-1S	2.5 cm x 5 m skin	box of 1
25288	1533NP-1SD	2.5 cm x 5 m skin with dispenser	box of 1
25290	1530NP-1S	2.5 cm x 5 m white Swiss design	box of 1
25291	1530NP-0SD	1.25 cm x 5 m white with dispenser	box of 1
34893	1530NP-1SD	h 2.5 cm x 5 m white with dispenser	box of 1

GIMA code	3M™ code	TRANSPORE™ - WHITE STERILE	Minimum order
25292	1527NP-1S	2.5 cm x 5 m	box of 1
25293	1527NP-1SD	2.5 cm x 5 m with dispenser	box of 1


OPTICLUDE™ - EYE PATCHES

Coloured eye patches for orthoptic therapies made of non-woven fabric, available in mini or maxi sizes for boys and girls.

Easy to apply and remove, without causing pain or discomfort.

Breathable, hypoallergenic and latex free.

GIMA code	3M™ code	OPTICLUDE™ - STERILE	Minimum order
25296	2537PE	OPTICLUDE™ eye patch mini - 5x6 cm	box of 30
25297	2539PE	OPTICLUDE™ eye patch maxi - 5.7x8 cm	box of 30



25296

Multilingual instructions: GB, FR, IT, ES, PT, DE, NL, DK, SE, NO, FI, LV, EE, PL, CZ, SK, HU, SI, HR, RO, BG, GR, TR.



25297