

WUNDER BODY COMPOSITION ANALYSER WITH PRINTER - Class III

25015 WUNDER WBA300 BODY COMPOSITION ANALYSER WITH PRINTER - 300 kg - Class III

Professional electronic scale with BIA body composition analyser. Simple, fast and non-invasive examination of BMI, fat, muscle and water components. Functions: weight, tare, mother/child weight, hold, BMI, autozero, programmable auto-off. - analysis of nutritional solutions advised by professionals - complete check-up in 20 sec. with a single weight reading - easy interpretation of results gives a complete picture, helping the patients to understand the importance of this screening. Large painted aluminium base with 4 stainless steel tactile electrodes. ABS indicator with LCD display set on a column in painted aluminium. Supplied with power adaptor, RS232 interface, level, 4x large anti-slip adjustable feet, 2 integrated wheels for easy transport, thermal printer with its power cable, connector, 1 paper roll 57 x Ø 30 mm. Multilingual manual: GB, IT, ES.

- 25005 HEIGHT ROD 60-210 cm - optional
- 25017 PAPER ROLL 57 mm x 9m - box of 25 - spare



LCD Display with 5 digits 20mm (Weight), Multi-line LCD Display 3 digits 12mm (BMI / Age-Fat), Keyboard with 20 keys

MADE IN ITALY

Printer Data:
Date, Time, Body Type, Sex, Age, Height, Weight, Clothes Weight, BMI, BMR, Fat (kg -%), FFM, Muscle Mass, Protein Mass, Skeletal Muscle, Bone mineral content, TBW, Intra-extracellular water.

25015

TECHNICAL SPECIFICATIONS

Frequencies: 4 bipolar mono-frequency 50 kHz electrodes
Impedance frequency: 50kHz 500uA
Capacity: 300 kg
Age Range: 10 ~ 80 years, increase of 1 year
Size: 450x550xh 880 mm
Platform size: 450x340xh 85 mm
Weight: 11 kg
Power supply: adapter AC 12V 2A 100-240V 50/60 Hz recharg. battery
Power supply optional: 6x1.5V AAA batteries

Range: 200 ~ 10000 / 0.1g
Sensitivity: 100g/0.1%

Column size: 100x50xh 750 mm

25015

WUNDER PROFESSIONAL DIGITAL SCALES - CLASS III

- 25000 WUNDER RB200 DIGITAL SCALE - 250 kg - Class III - standard
- 25001 WUNDER RB200 DIGITAL SCALE - 250 kg - Class III - customized*
- 25003 WUNDER RA300 DIGITAL SCALE - 300 kg - Class III - standard

Professional digital column scales with a large and stable painted aluminium base with non-slip painted metal platform. The scale is equipped with a dual LCD display with triple simultaneous readout: weight, height and BMI. Adjustable display for clearer and more comfortable reading. Functions: weight, tare, mother/child weight, hold, BMI, autozero, programmable auto-off. 4 large adjustable feet and 2 integrated wheels for easy transport. Supplied with power adaptor, RS232 interface, level, 4 adjustable feet, 2 wheels and multilingual manual: GB, IT, ES. Available on gimaitaly.com: FR, PT, DE, NL, SE, NO, FI. MDR certified.

- 25005 HEIGHT ROD 60-210 cm - optional

*CUSTOMIZABLE OPTIONS for 25001 (available in 20 days)

- unit of measure kg/lb
- customized calibration based on gravity factor
- France: customized calibration by area with carnet métrologique

TECHNICAL SPECIFICATIONS

	25000-01	25003
Capacity:	150/250 kg	300 kg
Sensitivity:	50/100 g	100 g
Display:	double LCD: LCD 20mm weight with 5 active digits; LCD 12 mm height and BMI	
Scale size - weight:	275x430xh 900mm - 7 kg	355x570xh 950 mm - 9 kg
Packing size:	330x1100x290 mm	330x950x150 mm
Platform size:	275x295xh 55mm	355x310xh 75 mm
Column size:	Ø60 x h 800mm	Ø60xh 800 mm
Power supply (included):	adapter 12V - 100-240V, 50/60 HZ	
Power supply (optional):	rechargeable battery or 6xAAA alkaline batteries	

Double LCD display

25000/1/3 + 25005

CE M 0474

CE 0425

Class III

25003

Large, anti-slip platform

MADE IN ITALY

WUNDER PROFESSIONAL ELECTRONIC BABY SCALES

- 24990 WUNDER BABY02-1 PROFESSIONAL ELECTRONIC SCALE 6-15 kg - Class III
- 24991 WUNDER BABY02-2 COMPACT PROFESSIONAL ELECTRONIC SCALE 10-20 kg - Class III

Professional and portable electronic scale for babies, made of white ABS with digital indicator. Functions: weight, tare, milk intake function, hold, auto zero, programmable auto-off. Breast-milk-intake function determines how much milk the baby has suckled. Comfortable and safe surface thanks to the curved edges. Easy to read double display.

MDR certified. Supplied with power supply, RS232 interface, level, adjustable feet and safety brackets.

- 24992 MANUAL HEIGHT ROD WITH ANALOGUE READING 35-80 cm - optional Height rod size: 87x30xh 7 cm

MADE IN ITALY

24991 + 24992

User manual: GB, IT, ES.
Available on gimaitaly.com:
FR, PT, DE, NL, SE, NO, FI.

CE 0425

CE M 0474

Class III

24990 + 24992

Telescopic, aluminium height rod with double ABS paddle (fixed headrest and sliding footrest)

Anti-overturning bars

TECHNICAL SPECIFICATIONS

	24990	24991
Capacity - sensitivity:	6-15 kg - 2/5 g (24990); 10-20 kg - 5/10 g (24991)	560x450xh 460 mm (24990); 560x540xh 145 mm (24991) - 4,8 kg
Display:	LCD 20mm weight with 5 active digits	adapter 12V - 2 A
Weighing surface:	LCD 12mm milk intake	rechargeable battery or 6xAAA alkaline batteries
Column height:	560x290 mm	480x640xh 910 mm
	460 mm	Cart load capacity:
		100 kg